




DENTAL INSTRUMENTS CATALOGUE



MAST PAK SURGICAL CORPORATION

Ugoki Road, Opposite DHL Main Office

Sialkot - 51310, Pakistan

Ph: +92 3248609456

Email: mirza.mastapak@gmail.com

Web: www.mastpaksurgicalcorp.com

TABLE OF CONTENT

01

FORCEPS

- » Extraction Forceps
- » Root Fragment Forceps

40

ELEVATORS

- » Root and Luxating Elevators

71

SCALPELS

- » Scalpel Handles and Blades

78

SCISSORS

- » Micro Scissors
- » Surgical Scissors
- » Ligature Scissors
- » Wire and Plate

106

TWEEZERS

- » Micro Forceps
- » Titanium Forceps
- » Dressing Forceps
- » Endo Forceps
- » Tissue Forceps
- » TC Forceps
- » Membrane

133

DRIVERS

- » Needle Holders

155

HEMOSTATS / GRASPING FORCEPS

- » Towel Clamps
- » Sponge and Dressing Forceps
- » Hemostatic Forceps
- » Tissue Grasping Forceps

166

RETRACTORS

- » Retractors
- » Mouth Gags and Tongue Depressors
- » Mouth Props

187

IMPLANTOLOGY

- » Oteotomes and Chisels
- » Mallets
- » Periosteals
- » Trephines
- » Periosteals
- » Osteotomes
- » Trephines
- » Bone and Tissue Harvest Accessories
- » Measuring Instruments
- » Bone Holding Forceps

210

PERIOSTEALS / CURETTES / BONE

- » Periosteal Elevators
- » Curettes
- » Bone Files
- » Tissue Nippers
- » Bone Rongeurs
- » Bone Cutting Forceps

235

PERIODONTICS / ENDODONTICS / RESTORATION

- » Periodontics
- » Endodontics
- » Diagnostic
- » Curettes and Scalers
- » Amalgam
- » Composite
- » Spatulas
- » Filling
- » Dental Restoration
- » Perio Chisels & Knives
- » Perio Probes & Explorers
- » Filling
- » Curettes and Scalers
- » Rubber-Dam

300

CROWN & BRIDGE

- » Crown and Bridge Instruments

Extracting Forceps, English Pattern

Zahnzangen, Englische Form

Pinze per estrazione, Forma Inglese

Pinzas da extracción, Forma Inglesa

Daviers, Forme Anglaise


Forceps

Pinzetten

Pinzas

Pinces


Pinzei



1

Fig. 1


Obere Schneide- und Eckzähne
Upper Incisors and cuspids
Incisivi e canini superiori
Incisivos y caninos superiores
Incisives et canines supérieures



2

Fig. 2


Obere Schneide- und Eckzähne
Upper Incisors and cuspids
Incisivi e canini superiori
Incisivos y caninos superiores
Incisives et canines supérieures



7

Fig. 7


Obere Prämolaren
Upper premolars
Premolari superiori
Premolares superiores
Prémolaires supérieures



17

Fig. 17

Obere Molaren, rechts
Upper molars, right
Molari superiori, destra
Molares superiores, derechos
Molaires supérieures, droites



18

Fig. 18

Obere Molaren, links
Upper molars, left
Molari superiori, sinistra
Molares superiores, izquierdos
Molaires supérieures, gauches



Extracting Forceps, English Pattern

Zahnzangen, Englische Form

Pinze per estrazione, Forma Inglese

Pinzas da extracci6n, Forma Inglesa



Daviers, Forme Anglaise



29

Fig. 29



Obere Wurzeln und Schneidezähne
Upper Roots and incisors
Radici e incisivi superiori
Raíces y incisivos superiores
Racines et incisives supérieures



30

Fig. 30

Obere Wurzeln
Upper Roots
Radici superiori
Raíces superiores
Racines supérieures



76S

Fig. 76

Obere Schneidezähne, Prämolaren, Wurzeln
Upper incisors, premolars, roots
Incisivi, premolari e radici superiori
Incisivos, premolares y raíces superiores
Incisives, premolaires et racines supérieures

Extracting Forceps, English Pattern

Zahnzangen, Englische Form

Pinze per estrazione, Forma Inglese

Pinzas da extracción, Forma Inglesa

Daviers, Forme Anglaise


Forceps

Pinzetten

Pinzas

Pinces


Pinzei



51

Fig. 51


Obere Wurzeln
Upper Roots
Radici superiori
Raíces superiores
Racines supérieures



51A

Fig. 51A


Obere Wurzeln
Upper Roots
Radici superiori
Raíces superiores
Racines supérieures



51L

Fig. 51AL

Obere Wurzeln
Upper Roots
Radici superiori
Raíces superiores
Racines supérieures



67A

Fig. 67A

Obere Dritte Molaren - Für Obere Weisheitszähne
Upper molars - For Upper wisdom teeth
Ottavi superiori - Per denti del giudizio superiori
Molares superiores - Para muelas del juicio superiores
Molaires supérieures - Pour dents de sagesse supérieures


Extracting Forceps, English Pattern

Zahnzangen, Englische Form

Pinze per estrazione, Forma Inglese


Pinzas da extracción, Forma Inglesa

Daviers, Forme Anglaise




33
Fig. 33

Untere Prämolaren, Schneidezähne und Wurzeln
Lower incisors, premolars and roots
Incisivi, premolari e radici inferiori
Incisivos, premolares y raíces inferiores
Incisives, prémolaires et racines inférieures



33A
Fig. 33A


Untere Wurzeln, schlankes Maul
Lower roots, Slender jaw
Radici inferiori, morso stretto
Raíces inferiores, boca estrecha
Racines inférieures, avec un bec très fin



13


Fig. 13

Untere Eckzähne und Prämolaren
Lower cuspids and premolars
Canini e premolari inferiori
Caninos y premolares inferiores
Canines et prémolaires inférieures



22
Fig. 22

Untere Molaren
Lower molars
Molari inferiori
Molares inferiores
Molaires inférieures



73
Fig. 73

Untere Molaren
Lower molars
Molari inferiori
Molares inferiores
Molaires inférieures

Extracting Forceps, English Pattern

Zahnzangen, Englische Form

Pinze per estrazione, Forma Inglese

Pinzas da extracción, Forma Inglesa

Daviers, Forme Anglaise





Forceps

Pinzetten

Pinzas



Pinces

Pinzei





74N
Fig. 74N

Untere Wurzeln, schneide-und Eckzähne
Lower roots, Incisors and bicuspid
Radici, incisivi e canini inferiori
Raíces, incisivos y caninos inferiores
Racines, incisives et canines inférieures



75
Fig. 75

Untere Eckzähne und Prämolaren
Lower Cuspids and bicuspid
Canini e premolari inferiori
Caninos y premolares inferiores
Canines et prémolaires inférieures



86A
Fig. 86A

Untere Molaren
Lower molars
Molari inferiori
Molares inferiores
Molaires inférieures

74
Fig. 74
Untere Wurzeln
Lower roots
Radici inferiori
Raíces inferiores
Racines inférieures


Extracting Forceps, English Pattern

Zahnzangen, Englische Form

Pinze per estrazione, Forma Inglese

Pinzas da extracción, Forma Inglesa

Daviers, Forme Anglaise



4

Fig. 4


Untere Schneide- und Eckzähne

Lower Incisors and cuspids

Incisivi e canini inferiori

Incisivos y caninos inferiores

Incisives et canines inférieures



8

Fig. 8


Untere Prämolaren

Lower premolars

Premolari inferiori

Premolares inferiores

Prémolaires inférieures



21

Fig. 21


Untere Weisheitszähne und Molaren

Lower wisdom teeth and molars

Denti del giudizio e molari inferiori

Muelas del juicio y Molares inferiores

Dents de sagesse et molaires inférieures



79

Fig. 79


Untere Weisheitszähne

Lower wisdom teeth

Denti del giudizio inferiori

Muelas del juicio inferiores

Dent de sagesse inférieures



79A

Fig. 79A

Untere Molaren und Weisheitszähne

Lower molars and wisdom teeth

Molari e denti del giudizio inferiori

Molares y muelas del juicio inferiores

Molaires et dent de sagesse inférieures

Extracting Forceps, English Pattern

Zahnzangen, Englische Form

Pinze per estrazione, Forma Inglese

Pinzas da extracci6n, Forma Inglesa

Daviers, Forme Anglaise


Forceps

Pinzetten


Pinzas

Pinces


Pinzei




87
Fig. 87
Untere Molaren
Lower molars
Molari inferiori
Molares inferiores
Molaires inférieures



89
Fig. 89
Obere Molaren, rechts
Upper molars, right
Molari superiori, destra
Molares superiores, derechos
Molaires supérieures, droites



90
Fig. 90
Obere Molaren, links
Upper molars, left
Molari superiori, sinistra
Molares superiores, izquierdos
Molaires supérieures, gauches



AM-451
Fig. 451
Deep Grip

Untere Wurzeln, schneide-und Eckzähne
Lower roots, incisors and bicusps
Radici, incisivi e canini inferiori
Raíces, incisivos y caninos inferiores
Racines, incisives et canines inférieures


Extracting Forceps, English Pattern

Zahnzangen, Englische Form

Pinze per estrazione, Forma Inglese

Pinzas da extracci6n, Forma Inglesa


Daviers, Forme Anglaise



68

Fig. 68


Untere Wurzeln
Lower Roots
Radici inferiori
Raíces inferiores
Racines inférieures



68A

Fig. 68A


Untere Wurzeln
Lower Roots
Radici inferiori
Raíces inferiores
Racines inférieures



55

Fig. 55

Zum Trennen oberer Molarenwurzeln
For separating upper molar roots
Per radici dei molari superiori, presa in profondità
Para partir raíces de molares superiores, profundo
Pour séparer les racines des molaires supérieures



56

Fig. 56

Zum Trennen Unterer Molarenwurzeln
For separating lower molar roots
Per radici dei molari inferiori, presa in profondità
Para partir raíces de molares inferiores, profundo
Pour séparer les racines des molaires inférieures

Extracting Forceps, English Pattern

Zahnzangen, Englische Form

Pinze per estrazione, Forma Inglese

Pinzas da extracción, Forma Inglesa

Daviers, Forme Anglaise


Forceps

Pinzette

Pinzas


Pinces

Pinzei



ROUTURIER
22 ½ R

Untere Molaren und Weisheitszähne, rechts
Lower molars and wisdom teeth, right
Molari e denti del giudizio inferiori, a destra
Molares y muelas del juicio inferiores, derecha
Molaires et dent de sagesse inférieures, droites



ROUTURIER
22 ½ L

Untere Molaren und Weisheitszähne, links
Lower molars and wisdom teeth, left
Molari e denti del giudizio inferiori, a sinistra
Molares y muelas del juicio inferiores, izquierda
Molaires et dent de sagesse inférieures, gauches


Extracting Forceps, English Pattern

Zahnzangen, Englische Form

Pinze per estrazione, Forma Inglese


Pinzas de extracción, Forma Inglesa

Daviers, Forme Anglaise




49
Fig. 49

Obere Wurzeln sehr feine
Upper very fine roots
Radici superiori molto sottili
Raíces superiores muy finas
Racines supérieures très fines



48
Fig. 48 - 44

Obere Wurzeln
Upper roots
Radici superiori
Raíces superiores
Racines supérieures



46L
Fig. 46L

Untere Wurzeln
Lower roots
Radici inferiori
Raíces inferiores
Racines inférieures

Extracting Forceps, Anatomical Handle

Zahnzangen, Paßformgriff

Pinze per estrazione, Manico anatomico

Pinzas da extracción, Mango anatómico

Daviers, Manche anatomique


Forceps

Pinzetten

Pinzas


Pinces

Pinzei




S-3AN
Fig. 3

Obere Schneide- und Eckzähne
Upper incisors and cuspids
Incisivi e canini superiori
Incisivos y caninos superiores
Incisives et canines supérieures




S-34
Fig. 34

Obere Schneide- und Eckzähne
Upper incisors and cuspids
Incisivi e canini superiori
Incisivos y caninos superiores
Incisives et canines supérieures



S-7
Fig. 7

Obere Prämolaren
Upper premolars
Premolari superiori
Premolares superiores
Prémolaires supérieures



S-35
Fig. 35

Obere Schneidezähne und Prämolaren, tiefgreifend
Upper incisors and bicuspid, deep-gripping
Incisivi e premolari superiori, presa in profondità
Incisivos y premolares superiores, a gran calado
Incisives et prémolaires supérieures, à attaque profonde


Extracting Forceps, Anatomical Handle

Zahnzangen, Paßformgriff


Pinze per estrazione, Manico anatomico

Pinzas da extracción, Mango anatómico


Daviers, Manche anatomique




S-31
Fig. 31



1/1




1/1




S-66R

Fig. 66R (89)


Obere Molaren, rechts
Upper molars, right
Molari superiori, destra
Molares superiores, derechos
Molaires supérieures, droites



1/1




1/1



S-66L

Fig. 66L (90)

Obere Molaren, links
Upper molars, left
Molari superiori, sinistra
Molares superiores, izquierdos
Molaires supérieures, gauches



1/1

Obere Prämolaren
Upper premolars
Premolari superiori
Premolares superiores
Prémolaires supérieures

Extracting Forceps, Anatomical Handle

Zahnzangen, Paßformgriff

Pinze per estrazione, Manico anatomico

Pinzas da extracción, Mango anatómico

Daviers, Manche anatomique


Forceps

Pinzetten

Pinzas


Pinces

Pinzei




S-67R
Fig. 67R

Obere Molaren, rechts
Upper molars, right
Molari superiori, destra
Molares superiores, derechos
Molaires supérieures, droites



S-67L
Fig. 67L

Obere Molaren, links
Upper molars, left
Molari superiori, sinistra
Molares superiores, izquierdos
Molaires supérieures, gauches



S-67N
Fig. 67N

Obere Molaren
Upper molars
Molari superiori
Molares superiores
Molaires supérieures


Extracting Forceps, Anatomical Handle

Zahnzangen, Paßformgriff


Pinze per estrazione, Manico anatomico

Pinzas da extracción, Mango anatómico


Daviers, Manche anatomique



FELSCH
S-80
Fig. 80




1/1




S-51B
Fig. 51B

Obere Dritte Molaren
Upper molars
Ottavi superiori
Molares superiores
Molaires supérieures



S-97
Fig. 97

Obere Wurzelplitterzange
Upper Root Fragments
Frammenti e radici superiori
Fragmentos de raíces superiores
Fragments des racines supérieures



S-97B
Fig. 97B

Obere Wurzelplitterzange
Upper Root Fragments
Frammenti e radici superiori
Fragmentos de raíces superiores
Fragments des racines supérieures

Für obere Weisheitszähne, mit Fangvorrichtung
For upper wisdom teeth, with retaining device
Denti del giudizio superiori, con dispositivo di presa
Muelas del juicio superiores, con dispositivo de retencion
Pour dents de sagesse supérieures, mors conçus pour recueillir la dent

Extracting Forceps, Anatomical Handle

Zahnzangen, Paßformgriff

Pinze per estrazione, Manico anatomico

Pinzas da extracción, Mango anatómico

Daviers, Manche anatomique

Forceps


Pinzetten

Pinzas


Pinces

Pinzei


Obere Wurzeln sehr feine
Upper very fine roots
Radici superiori molto sottili
Raíces superiores muy finas
Racines supérieures très fines



S-49
Fig. 49



S-51LX
Fig. 51LX



S-351S
Fig. 351S


Extracting Forceps, Anatomical Handle

Zahnzangen, Paßformgriff

Pinze per estrazione, Manico anatomico

Pinzas da extracción, Mango anatómico


Daviers, Manche anatomique




I-74N

Fig. 74N

Untere Wurzeln Schneide- und Eckzähne
Lower roots, incisors and cuspids
Radici, incisivi e canini inferiori
Raíces, incisivos y caninos inferiores
Racines, incisives et canines inférieures




I-74
Fig. 74



I-36

Fig. 36

Untere Wurzeln Schneide- und Eckzähne
Lower roots, incisors and cuspids
Radici, incisivi e canini inferiori
Raíces, incisivos y caninos inferiores
Racines, incisives et canines inférieures



I-36N

Fig. 36N

Untere Wurzeln Schneide- und Eckzähne
Lower roots, incisors and cuspids
Radici, incisivi e canini inferiori
Raíces, incisivos y caninos inferiores
Racines, incisives et canines inférieures

Extracting Forceps, Anatomical Handle

Zahnzangen, Paßformgriff

Pinze per estrazione, Manico anatomico

Pinzas da extracción, Mango anatómico

Daviers, Manche anatomique

Forceps

Pinzette

Pinzas

Pinces

Pinzei




I-59
Fig. 59

Untere Wurzelplitterzange
Lower Root Fragments
Frammenti e radici inferiori
Fragmentos de raíces inferiores
Fragments des racines inférieures




I-86C
Fig. 86C

Untere Molaren mit kariesen oder abgebrochenen Kronen
Lower molars with carious or broken caps
Molari inferiori con corone cariate o spezzate
Molares inferiores con coronas rotas o cariadas
Molaires inférieures avec couronnes cariées ou cassées




I-73
Fig. 73

1/1




I-13
Fig. 13

Untere Prämolaren und Eckzähne
Lower premolars and cuspids
Canini e premolari inferiori
Caninos y premolares inferiores
Canines et prémolaires inférieures



1/1



I-22
Fig. 22

Untere Molaren
Lower molars
Molari inferiori
Molares inferiores
Molaires inférieures



1/1

Untere Molaren
Lower molars
Molari inferiori
Molares inferiores
Molaires inférieures

Extracting Forceps, Anatomical Handle

Zahnzangen, Paßformgriff

Pinze per estrazione, Manico anatomico

Pinzas da extracción, Mango anatómico


Daviers, Manche anatomique




Untere Weisheitszähne
Lower wisdom teeth
Denti del giudizio inferiori
Muelas del juicio inferiores
Dent de sagesse inférieures

Untere Prämolaren
Lower premolars
Premolari inferiori
Premolares inferiores
Prémolaires inférieures


Untere Molaren
Lower Molars
Molari inferiori
Molares inferiores
Molaires inférieures




I-79
Fig. 79



I-79N
Fig. 79N



I-46N
Fig. 46N



I-87
Fig. 87

Extracting Forceps, Anatomical Handle

Zahnzangen, Paßformgriff

Pinze per estrazione, Manico anatomico

Pinzas da extracción, Mango anatómico

Daviers, Manche anatomique

Forceps

Pinzetten


Pinzas

Pinces

Pinzei



I-46LX
Fig. 46LX



I-46L
Fig. 46L

Extracting Forceps, Anatomical Handle

Zahnzangen, Paßformgriff



Pinze per estrazione, Manico anatomico

Pinzas da extracción, Mango anatómico

Daviers, Manche anatomique




Untere Molaren und Weisheitszähne
Lower molars and wisdom teeth
Molari e denti del giudizio inferiori
Molares y muelas del juicio inferiores
Molaires et dent de sagesse inférieures



Links
Left
Sinistra
Izquierdos
Gauches

I-22 1/2 L



Rechts
Right
Destra
Derechos
Droites

I-22 1/2 R

Extracting Forceps, Anatomical Handle

Zahnzangen, Paßformgriff

Pinze per estrazione, Manico anatomico

Pinzas da extracción, Mango anatómico

Daviers, Manche anatomique

Forceps


Pinzetten

Pinzas

Pinces


Pinzei

Untere Wurzeln sehr feine
Lower very fine roots
Radici inferiori molto sottili
Raíces inferiores muy finas
Racines inférieures très fines




Links
Left
Sinistra
Izquierdos
Gauches

I-133 L
Fig.133



Rechts
Right
Destra
Derechos
Droites


I-133 R
Fig.133



I-333S
Fig.333S


Extracting Forceps, American Pattern

Zahnzangen, Amerikanische Form
Pinze per estrazione, Forma Americana
Pinzas da extracción, Forma Americana
Daviers, Forme Américaine




AM-1
Fig. 1

Obere Schneide- und Eckzähne
Upper incisors and cuspids
Incisivi e canini superiori
Incisivos y caninos superiores
Incisives et canines supérieures




CRYER
AM-150
Fig. 150

Obere Prämolaren, Schneidezähne, Wurzeln
Upper Bicuspid, incisors and roots
Premolari, incisivi e radici superiori
Premolares, incisivos y raíces superiores
Prémolaires, incisives et racines supérieures



CRYER
AM-150A
Fig. 150A

Obere Prämolaren, Schneidezähne, Wurzeln
Upper Bicuspid, incisors and roots
Premolari, incisivi e radici superiori
Premolares, incisivos y raíces superiores
Prémolaires, incisives et racines supérieures



AM-62
Fig. 62



Schneidezähne, Prämolaren, Milchzähne
Incisors, premolars, deciduous teeth
Incisivi, premolari, denti da latte
Incisivos, premolares, dientes de leche
Incisives, prémolaires, dents de lait

Extracting Forceps, American Pattern

Zahnzangen, Amerikanische Form
Pinze per estrazione, Forma Americana
Pinzas da extracción, Forma Americana
Daviers, Forme Américaine



Forceps

Pinzetten
Pinzas
Pinces
Pinzei





AM-53R
Fig. 53R

Obere Molaren, Rechts
Upper Molars, Right
Molari superiori, Destra
Molares superiores, Derechos
Molaires supérieures, Droites





AM-53L
Fig. 53L

Obere Molaren, Links
Upper Molars, Left
Molari superiori, Sinistra
Molares superiores, Izquierdos
Molaires supérieures, Gauches



NEVIUS
AM-88R
Fig. 88R

Obere Molaren, Rechts
Upper Molars, Right
Molari superiori, Destra
Molares superiores, Derechos
Molaires supérieures, Droites




NEVIUS
AM-88L
Fig. 88L

Obere Molaren, Links
Upper Molars, Left
Molari superiori, Sinistra
Molares superiores, Izquierdos
Molaires supérieures, Gauches


Extracting Forceps, American Pattern

Zahnzangen, Amerikanische Form
Pinze per estrazione, Forma Americana
Pinzas da extracción, Forma Americana
Daviers, Forme Américaine




AM-65
Fig.65

Obere Schneidezähne und Wurzeln
Upper Incisors and Roots
Incisivi e radici superiori
Incisivos y raíces superiores
Incisives et racines supérieures



PARMLY
AM-32A
Fig. 32A

Obere Wurzeln und Prämolaren
Upper Roots and Premolars
Radici e premolari superiori
Raíces y premolares superiores
Racines et prémolaires supérieures



PARMLY
AM-32
Fig. 32


Obere Prämolaren und Molaren
Upper Premolars and Molars
Premolari e molari superiori
Premolares y Molares superiores
Prémolaires et Molaires supérieures

Extracting Forceps, American Pattern

Zahnzangen, Amerikanische Form
Pinze per estrazione, Forma Americana
Pinzas da extracción, Forma Americana
Daviers, Forme Américaine


Forceps

Pinzetten
Pinzas
Pinces
Pinzei




AM-24
Fig. 24


Obere Molaren
Upper Molars
Molari superiori
Molares superiores
Molaires supérieures



1/1 1/1




1/1 1/1



HARRIS
AM-18R
Fig. 18R

Obere Molaren, Rechts
Upper Molars, Right
Molari superiori, Destra
Molares superiores, Derechos
Molaires supérieures, Droites




1/1 1/1




HARRIS
AM-18L
Fig. 18L

Obere Molaren, Links
Upper Molars, Left
Molari superiori, Sinistra
Molares superiores, Isquierdos
Molaires supérieures, Gauches



1/1 1/1



AM-210
Fig. 210

Obere Molaren
Upper Molars
Molari superiori
Molares superiores
Molaires supérieures




Extracting Forceps, American Pattern

Zahnzangen, Amerikanische Form

Pinze per estrazione, Forma Americana

Pinzas da extracción, Forma Americana


Daviers, Forme Américaine



**TOMES
AM-69**
Fig.69

Obere und Untere
Upper and lower
Superiori e inferiori
Superiores y inferiores
Supérieures et inférieures


Wurzeln
Roots
Radici
Raíces
Racines



**CRYER
AM-151**
Fig.151

Untere
Lower
Inferiori
Inferiores
Inférieures

Schneidezähne, Prämolaren, Wurzeln
Incisors, Premolars, Roots
Incisivi, premolari, radici
Incisivos, premolares, raíces
Incisives, prémolaires, racines



**HULL
AM-101**
Fig. 101



Für Eckzähne und Prämolaren,
Universalzange und Kinderzange
For canines and Premolars,
Universal and Tooth forceps for children
Canini e premolari,
universale e per bambini
Caninos y premolares,
universal y para niños
Canines et prémolaires,
universel et pour enfants

Extracting Forceps, American Pattern

Zahnzangen, Amerikanische Form
Pinze per estrazione, Forma Americana
Pinzas da extracción, Forma Americana
Daviers, Forme Américaine



Forceps

Pinzetten
Pinzas
Pinces
Pinzei





AM-203
Fig. 203

Untere Schneidezähne, Prämolaren, Wurzeln
Lower Incisors, Premolars, Roots
Incisivi, premolari, radici inferiori
Incisivos, premolares, raíces inferiores
Incisives, prémolaires, racines inférieures



AM-17
Fig. 17

Untere Molaren
Lower Molars
Molari inferiori
Molares inferiores
Molaires inférieures



AM-23
Fig. 23

Untere Molaren
Lower Molars
Molari inferiori
Molares inferiores
Molaires inférieures


Extracting Forceps, American Pattern

Zahnzangen, Amerikanische Form

Pinze per estrazione, Forma Americana


Pinzas da extracción, Forma Americana

Daviers, Forme Américaine




AM-16
Fig. 16

Untere Molaren
Lower Molars
Molari inferiori
Molares inferiores
Molaires inférieures



HARRIS
AM-15
Fig. 15

Untere Molaren
Lower Molars
Molari inferiori
Molares inferiores
Molaires inférieures



AM-103
Fig. 103


Untere Schneidezähne, Prämolaren, Wurzeln
Lower Incisors, Premolars, Roots
Incisivi, premolari, radici inferiori
Incisivos, premolares, raíces inferiores
Incisives, prémolaires, racines inférieures

Extracting Forceps, American Pattern

Zahnzangen, Amerikanische Form
Pinze per estrazione, Forma Americana
Pinzas da extracción, Forma Americana
Daviers, Forme Américaine


Forceps

Pinzetten
Pinzas
Pinces
Pinzei




AM-222
Fig. 222

Untere Weisheitszähne
Lower wisdom teeth
Denti del giudizio inferiori
Muelas del juicio inferiores
Dent de sagesse inférieures



COHEN
AM-250

Universal
Universal
Universale
Universal
Universeel






PHYSICK
AM-5
Fig. 5

Untere Weisheitszähne
Lower wisdom teeth
Denti del giudizio inferiori
Muelas del juicio inferiores
Dent de sagesse inférieures

Extracting Forceps, American Pattern



Zahnzangen, Amerikanische Form
Pinze per estrazione, Forma Americana
Pinzas da extracción, Forma Americana
Daviers, Forme Américaine



Obere
Upper
Superiori
Superiores
Supérieures

Prämolaren, Schneidezähne und Wurzeln
Bicuspids, incisors and roots
Premolari, incisivi e radici
Premolares, incisivos y raíces
Premolaires, incisives et racines

AM-150X
Fig. 150AS



Untere
Lower
Inferiori
Inferiores
Inférieures

Prämolaren, Schneidezähne und Wurzeln
Bicuspids, incisors and roots
Premolari, incisivi e radici
Premolares, incisivos y raíces
Premolaires, incisives et racines

AM-151X
Fig. 151AS


Extracting Forceps, American Pattern

Zahnzangen, Amerikanische Form
Pinze per estrazione, Forma Americana
Pinzas da extracción, Forma Americana
Daviers, Forme Américaine

Forceps

Pinzetten
Pinzas
Pinces
Pinzei

Für Kinder
For children
Per bambini
Para niños
Pour enfants



Extracting Forceps, American Pattern

Zahnzangen, Amerikanische Form

Pinze per estrazione, Forma Americana

Pinzas da extracción, Forma Americana

Daviers, Forme Américaine



Geriefte - Serrated - Striata - Estriadas - Striées


Obere Schneidezähne, Prämolaren, Wurzeln

Upper incisors, premolars, roots

Incisivi, premolari e radici superiori

Incisivos, premolares y raíces superiores

Incisives, premolaires et racines supérieures





Untere Schneidezähne, Prämolaren, Wurzeln

Lower incisors, premolars, roots


Incisivi, premolari e radici inferiori

Incisivos, premolares y raíces inferiores

Incisives, premolaires et racines inférieures




CRYER
AMG-150S
Fig.150 S




CRYER
AMG-150
Fig. 150

Für Erwachsene
For adults
Per adulti
Para adultos
Pour adults

Für Kinder
For children
Per bambini
Para niños
Pour enfants



CRYER
AMG-151S
Fig.151 S



CRYER
AMG-151
Fig. 151

Für Erwachsene
For adults
Per adulti
Para adultos
Pour adults

Für Kinder
For children
Per bambini
Para niños
Pour enfants

Extracting Forceps

Zahnzangen

Pinze per estrazione

Pinzas da extracción

Daviers

Forceps

Pinzetten

Pinzas

Pinces

Pinzei

Schneidezähne, Eckzähne, Wurzeln

Incisors, bicuspid and roots

Incisivi, canini e radici

Incisivos, caninos y raíces

Incisives, canines et racines


Molaren

Molars

Molari


Molares

Molaires



MEAD
MD-1

MEAD
MD-3



MEAD
MD-2

MEAD
MD-4

Extracting Forceps for children, American Pattern

Kinderzahnzangen, Amerikanische Form

Pinze per estrazione per bambini, Forma Americana

Forceps dentales para niños, Forma Americana

Daviers pour enfants, Forme Américaine




Molaren
Pedo Lower Molars
Molari
Molares
Molaires


Molaren
Pedo Lower Molars
Molari
Molares
Molaires

Zähne und wurzeln
Pedo upper incisors
Bicuspids and roots
Denti e radici
Dientes y raíces
Dentes et racines


Zähne und wurzeln
Pedo lower incisors
Bicuspids and roots
Denti e radici
Dientes y raíces
Dentes et racines




PEB-17S
Fig. 17SK



PEB-23S
Fig. 23SK



PEB-150S
Fig. 150SK




PEB-151S
Fig. 151SK

Extracting Forceps for children

Kinderzahnzangen
Pinze per estrazione per bambini
Forceps dentales para niños
Daviers pour enfants


Forceps

Pinzetten
Pinzas
Pinces
Pinzei




PEB-1
Fig. 137

Obere Schneidezähne
Upper Incisors
Incisivi superiori
Incisivos superiores
Incisives supérieures




PEB-2
Fig. 139

Obere Prämolaren
Upper Premolars
Premolari superiori
Premolares superiores
Prémolaires supérieures




PEB-4
Fig. 51S

Obere Wurzeln
Upper roots
Radici superiori
Raíces superiores
Racines supérieures




PEB-5
Fig. 5

Untere Schneidezähne
Lower incisors
Incisivi inferiori
Incisivos inferiores
Incisives inférieures




PEB-3

Obere Molaren
Upper molars
Molari superiori
Molares superiores
Molaires supérieures



PEB-6
Fig. 6

Untere Molaren
Lower Molars
Molari inferiori
Molares inferiores
Molaires inférieures



PEB-7
Fig. 7

Untere Wurzeln
Lower roots
Radici inferiori
Raíces inferiores
Racines inférieures


Root Fragment Forceps

Wurzelfragmentinstrumente

Strumenti per frammenti di radici


Instrumentos para fragmentos de raíces

Instruments pour fragments de racines




Obere Wurzeln, Universal
Upper roots, Universal
Radici superiori, Universale
Raíces superiores, Universal
Racines superieures, Universel

**WITZEL
DK-26**
170 mm




Untere Wurzeln, Universal
Lower roots, Universal
Radici inferiori, Universale
Raíces inferiores, Universal
Racines inferieures, Universel

**WITZEL
DK-28**
165 mm



Obere Wurzeln, Universal
Upper roots, Universal
Radici superiori, Universale
Raíces superiores, Universal
Racines superieures, Universel

S-400
400 (300)
165 mm



Untere Wurzeln, Universal
Lower roots, Universal
Radici inferiori, Universale
Raíces inferiores, Universal
Racines inferieures, Universel

I-401
401 (301)
160 mm

Root Fragment Forceps

Wurzelfragmentinstrumente

Strumenti per frammenti di radici

Instrumentos para fragmentos de raíces

Instruments pour fragments de racines


Forceps

Pinzetten


Pinzas

Pinces

Pinzei




PEET
DR-40
120 mm




RO-1
175 mm

Untere Molaren und Weisheitszähne
Lower molars and wisdom teeth
Denti del giudizio e molari inferiori
Muelas del juicio y molares inferiores
Dent de sagesse et molaires inférieures

Wurzelfragmente
Root fragments
Frammenti di radici
Fragmentos de raíces
Fragments de racines



STIEGLITZ
DR-50/DR-60
140 mm



STIEGLITZ
DR-70
120 mm


Root Fragment Forceps

Wurzelfragmentinstrumente

Strumenti per frammenti di radici


Instrumentos para fragmentos de raíces

Instruments pour fragments de racines




DA-01M
130 mm

Diamantiert
With Diamond
Con diamante
Diamantada
Diamantées




DA-10
130 mm


Diamantiert
With Diamond
Con diamante
Diamantada
Diamantées



DA-01
130 mm



DA-02
130 mm



DA-03
125 mm

ELEVATORS




Root Elevators

Wurzelheber

Leve per radici


Elevadores para raíces

Elevateurs à racines



SMALL HANDLE

Luxierinstrumente | Luxating Elevators | Elevatori di lussazione | Elevadores de luxación | Élévateurs de luxation



Root Elevators








Wurzelheber
Leve per radici
Elevadores para raíces
Elevateurs à racines

Elevators

Aufzüge
Ascenseurs
Ascensori
Lifty



DL-81/DL-87
150 mm



Root Elevators

Wurzelheber


Leve per radici

Elevadores para raíces

Elevateurs à racines



2,0 mm




DL-91
150 mm

2,5 mm




DL-92
150 mm

3,0 mm




DL-93
150 mm

3,5 mm




DL-94
150 mm

4,0 mm




DL-95
150 mm

4,5 mm




DL-96
150 mm

5,0 mm



DL-97
150 mm



DL-91/DL-97
150 mm

Root Elevators





Wurzelheber
Leve per radici
Elevadores para raíces
Élevateurs à racines

Elevators




Aufzüge
Ascenseurs
Ascensori
Lifty

SMALL HANDLE



Luxierinstrumente | Luxating Elevators | Elevatori di lussazione | Elevadores de luxación | Élévateurs de luxation



DL-102/DL-105
150 mm



DL-112/DL-115
150 mm



Root Elevators

Wurzelheber

Leve per radici


Elevadores para raíces

Elevateurs à racines




Wurzelspitzenheber
Root fragment Elevators
Leve per frammenti di radice
Elevadores para fragmentos de raíces
Élevateurs de fragments de racines


HEIDBRINK
HE-2S
Right # 2
Root Tip Picks
140 mm
SMALL HANDLE





HEIDBRINK
HE-3S
Left # 3
Root Tip Picks
140 mm
SMALL HANDLE




Gezahnt
Serrated
Dentellata
Dentada
Dentelé





HEIDBRINK
HE-1
Straight # 1
Root Tip Picks
140 mm




HEIDBRINK
HE-2
Right # 2
Root Tip Picks
140 mm




HEIDBRINK
HE-3
Left # 3
Root Tip Picks
140 mm




HEIDBRINK
725-1/2
195 mm





HEIDBRINK
725-3/4
195 mm



HEIDBRINK
725-5/6
195 mm



HEIDBRINK
725-3/4S
195 mm



Root Elevators


Wurzelheber
Leve per radici
Elevadores para raíces
Élevateurs à racines

Elevators

Aufzüge
Ascenseurs
Ascensori
Lifty

Extraktionshebel für obere Weisheitszähne
Root Elevator for upper wisdom teeth
Leva per estrazione dente del giudizio superiore
Elevador de extracción de muelas del juicio superiores
Élévateur à racines pour les dents de sagesse supérieures

Wurzelheber
Root Elevators
Leve per radici
Elevadores para raíces
Élevateurs de racines



ROSIN
DL-58
160 mm



WARWICK-JAMES
DL-20
160 mm

Obere
Upper
Superiori
Superiores
Supérieures



WARWICK-JAMES
DL-21
155 mm

Untere
Lower
Inferiori
Inferiores
Inférieures



WARWICK-JAMES
DL-22
155 mm

Untere
Lower
Inferiori
Inferiores
Inférieures

Root Elevators

Wurzelheber

Leve per radici

Elevadores para raíces

Elevateurs à racines



Eine seite gezahnt
One side serrated
Un lato dentellato
Dentadura mäs lado
Dentelure sur un côté






Fig. 94



DL-76
155 mm

Fig. 95



DL-77
155 mm

Fig. 100L



DL-78
155 mm

Fig. 100R



DL-79
155 mm



Root Elevators

Wurzelheber
Leve per radici
Elevadores para raíces
Elevateurs à racines


Elevators

Aufzüge
Ascenseurs
Ascensori
Lifty


Gezahnt
Serrated
Dentellate
Dentadas
Dentelé




APEXO
DL-34S
301
155 mm




APEXO
DL-34S
301
130 mm
SMALL HANDLE




APEXO
DL-35/DL-36
155 mm




DL-35




DL-36



APEXO
DL-37/DL-38
155 mm



DL-37



DL-38


Root Elevators

Wurzelheber

Leve per radici

Elevadores para raíces

Elevateurs à racines



BEIN-1T
150 mm

BEIN-2
155 mm

BEIN-3
155 mm

BEIN-4
155 mm

Root Elevators

Wurzelheber

Leve per radici

Elevadores para raíces

Elevateurs à racines


Elevators

Aufzüge


Ascenseurs

Ascensori


Lifty




FLOHR-1
140 mm




FLOHR-2
150 mm



FLOHR-3
150 mm




APICAL-1
155 mm




APICAL-4
155 mm



1/1



APICAL-2
155 mm



APICAL-3
155 mm

Root Elevators

Wurzelheber


Leve per radici

Elevadores para raíces

Elevateurs à racines




Fig. 44 1/1




CRYER
DL-27
155 mm

Fig. 45 1/1




CRYER
DL-28
155 mm

Fig. 25 1/1




CRYER-DX
DL-32
155 mm


Fig. 26 1/1



CRYER-SX
DL-33
155 mm




CRYER
155 mm



KRALLENHEBER
DL-1
135 mm

Untere Wurzeln
Lower Roots
Radici inferiori
Racines inférieures
Raíces inferiores



KRALLENHEBER
DL-2
135 mm



Untere Wurzeln
Lower Roots
Radici inferiori
Racines inférieures
Raíces inferiores

Root Elevators

Wurzelheber
Leve per radici
Elevadores para raíces
Elevateurs à racines

Elevators

Aufzüge
Ascenseurs
Ascensori
Lifty




Root Elevators

Wurzelheber



Leve per radici

Elevadores para raíces


Elevateurs à racines



MILLER-APEXO
MI-10/MI-11
155 mm



MILLER-APEXO
MI-12
155 mm




MILLER-APEXO
MI-13
155 mm

Root Elevators

Wurzelheber
Leve per radici
Elevadores para raíces
Elevateurs à racines

Elevators

Aufzüge
Ascenseurs
Ascensori
Lifty




Root Elevators

Wurzelheber


Leve per radici

Elevadores para raíces


Elevateurs à racines




**KOPP
DL-30**
130 mm




**KOPP
DL-31**
130 mm



BEIN-1Z
130 mm



BEIN-2Z
155 mm




**LINDO-LEVIEN
BEIN-4Z**
155 mm

Root Elevators

Wurzelheber
Leve per radici
Elevadores para raíces
Elevateurs à racines

Elevators

Aufzüge
Ascenseurs
Ascensori
Lifty



Root Elevators

Wurzelheber


Leve per radici

Elevadores para raíces


Elevateurs à racines




SELDIN
DL-48
304W
155 mm



SELDIN
DL-49
34S
155 mm




304W





SELDIN
DL-48
304W
155 mm



34S



SELDIN
DL-49
34S
155 mm



DL-54
160 mm



DL-55
160 mm

Root Elevators

Wurzelheber

Leve per radici

Elevadores para raíces

Elevateurs à racines

Elevators


Aufzüge

Ascenseurs

Ascensori


Lifty

Fig. 6




POTT
DL-98
160 mm

Fig. 7



POTT
DL-99
160 mm



POTT
DL-24/DL-25
155 mm

Root Elevators

Wurzelheber


Leve per radici

Elevadores para raíces

Elevateurs à racines




Syndesmotome
Syndesmotomes
Sindesmotomi
Syndesmôtomos
Syndesmotomes




BERNARD
DL-73/DL-124

Gerade
Straight
Diritto
Rectas
Droites




BERNARD
DL-73
150 mm

Links
Left
Sinistra
A la izquierda
À gauche




BERNARD
DL-74
140 mm

Rechts
Right
Destro
A la derecha
À droite




BERNARD
DL-75
140 mm

Bajonett - Bayonet
Baionetta - Bayoneta
Baionnette



BERNARD
DL-123
150 mm




BERNARD
DL-124
150 mm

Root Elevators

Wurzelheber
Leve per radici
Elevadores para raíces
Elevateurs à racines

Elevators


Aufzüge
Ascenseurs
Ascensori
Lifty



Meissel
Chisel
Scalpello
Cincele
Ciseaux-burins

Gerade
Straight
Diritto
Recto
Droit

DL-88
130 mm



Meissel
Chisel
Scalpello
Cincele
Ciseaux-burins


Bajonett
Bayonet
Baionetta
Bayoneta
Baionnette

DL-89
130 mm


Syndesmotome
Syndesmotomes
Sindesmotomo
Syndesmòtomos
Syndesmotomes



DL-63



DL-62




DL-61


Raspatorien
Raspatories
Scollaperiosteo
Legras
Rugines



DL-59



DL-68



DL-69


Root Elevators

Wurzelheber

Leve per radici

Elevadores para raíces


Elevateurs à racines



Untere Wurzeln
Lower Roots
Radici inferiori
Raíces inferiores
Racines inférieures

Links
Left
Sinistra
Izquierda
Gauche

WINTER
DL-100
Fig.1



Untere Wurzeln
Lower Roots
Radici inferiori
Raíces inferiores
Racines inférieures

Rechts
Right
Dextra
Derecha
Droites

WINTER
DL-101
Fig.2

Root Elevators

Wurzelheber
Leve per radici
Elevadores para raíces
Elevateurs à racines

Elevators


Aufzüge
Ascenseurs
Ascensori
Lifty

Fig. 1L



WINTER
DL-64
120 mm

Fig. 1R



WINTER
DL-65
120 mm

Fig. 14L



WINTER
DL-66
125 mm

Fig. 14R



WINTER
DL-67
125 mm


Root Elevators

Wurzelheber


Leve per radici

Elevadores para raíces

Elevateurs à racines




LECLUSE
DL-90
130 mm




BARRY
DL-110
Fig. 322

Links
Left
Sinistra
Izquierda
Gauche




BARRY
DL-108 Fig. 320
115 mm



BARRY
DL-111
Fig. 323

Rechts
Right
Destra
Derecha
Droites



BARRY
DL-109 Fig. 321
115 mm

Root Elevators

Wurzelheber
Leve per radici
Elevadores para raíces
Elevateurs à racines

Elevators

Aufzüge
Ascenseurs
Ascensori
Lifty

3 mm



LA-1
160 mm

3 mm



LA-2
160 mm

4 mm



LA-3
160 mm

4 mm



LA-4
160 mm

4,5 mm



LA-5
160 mm

4,5 mm



LA-6
160 mm

Root Elevators

Wurzelheber

Leve per radici

Elevadores para raíces

Elevateurs à racines



LA-7
160 mm



LA-8
160 mm



LA-9
160 mm



LA-10
160 mm

Root Elevators

Wurzelheber
Leve per radici
Elevadores para raíces
Elevateurs à racines

Elevators

Aufzüge
Ascenseurs
Ascensori
Lifty

2 mm



LA-11
155 mm

3 mm



LA-12
155 mm

4 mm



LA-13
155 mm



LA-16
160 mm



LA-17
160 mm

Root Elevators

Wurzelheber

Leve per radici

Elevadores para raíces

Elevateurs à racines



LA-20
160 mm



LA-21
160 mm



LA-22
160 mm



LA-23
160 mm

Syndesmotomes

Syndesmotome
Sindesmotomo
Syndesmòtomos
Syndesmotomes

Elevators

Aufzüge
Ascenseurs
Ascensori
Lifty



Massiv
Solid
Manico
Pieno
Massif

CHOMPRET
CH-11
150 mm



Massiv
Solid
Manico
Pieno
Massif

CHOMPRET
CH-12
150 mm



Massiv
Solid
Manico
Pieno
Massif

CHOMPRET
CH-13
150 mm



Massiv
Solid
Manico
Pieno
Massif

CHOMPRET
CH-14
145 mm



CHOMPRET
CH-1
165 mm



CHOMPRET
CH-2
165 mm



CHOMPRET
CH-3
165 mm



CHOMPRET
CH-4
160 mm

Scalers - Periotomes

Zahnreiniger - Periotome


Ablatori di tartaro - Periotomi

Limpiadores de sarro - Periotomos


Detarteurs - Périotomes



MORSE - SCALERS




Periotome - Periotomes
Periotomi - Periotomos
Périotomes



Ersatz
Insert for replacement
Inserto di ricambio
Puntas de recambio
Pointes de rechange

Periotome - Periotomes
Periotomi - Periotomos
Périotomes




Für / For
Per / Para / Pour
MA-9



MA-9R
48 mm

MA-9
Handle
140 mm
Rostfrei
Stainless Steel
Inox

MC-70
(MA-9 + L-71, L-72, L73)
Komplett mit 1 Griff + 3 Spatel
Complete with 1 Handle + 3 Spatulas
Completo di 1 Manico +3 Spatole
Completo con 1 mango + 3 Espátulas
Complete avec 1 Manche + 3 Spatules



Syndesmotomes

Syndesmotome
Sindesmotomo
Syndesmòtomos
Syndesmotomes

Elevators

Aufzüge
Ascenseurs
Ascensori
Lifty

SDM-C

Komplett mit 2 Griffe + 20 Klingen
Complete with 2 Handles + 20 Blades
Completo di 2 Manici + 20 lame
Completo con 2 mangos + 20 hojas
Complete avec 2 Manches + 20 lames

SD-M


Griff
Handle
Manico
Mango
Manche

SDM-9


Komplett mit 2 Griffe + 9 Klingen APICALI
Complete with 2 Handles + 9 Blades APICALI
Completo di 2 Manici + 9 lame APICALI
Completo con 2 mangos + 9 hojas APICALI
Complete avec 2 Manches + 9 lames APICALI

SDM-11

Komplett mit 2 Griffe + 11 Klingen CERVICALI
Complete with 2 Handles + 11 Blades CERVICALI
Completo di 2 Manici + 11 lame CERVICALI
Completo con 2 mangos + 11 hojas CERVICALI
Complete avec 2 Manches + 11 lames CERVICALI



SCALPELS



Scalpel Handles - Micro Scalpel Blades

Skalpellgriffe - Mikro Skalpellklingen

Manici per bisturi - Micro lame per bisturi

Mangos de bisturí - Micro hojas de bisturí

Manches pour bistouris - Micro lames de bistouri

Scalpels

Skalpelee

Bisturies

Bistouris

Bisturi





LB-40

Packung 10 Klingen
Package 10 blades
Scatola da 10 lame
Cajita de 10 hojas
Boite de 10 lames

Skalpellgriffe für Mikro-Klingen
Skalpel handles for Micro blades
Manici per micro lame
Mangos de bisturí para micro hojas
Manches de bistouris pour micro lames

Micro-Klingen
Micro blades
Micro lame
Micro hojas
Micro lames

Rasierklingen-Brecker und halter
Razor blade breaker and holder
Manico porta lame
Mango para Hojas de afeitar
Porte lames à rasoir



NK-12
120 mm



MA-11
135 mm



BB-46
135 mm



BB-48
155 mm



LB-361
Conf. 10 pz.



LB-362
Conf. 10 pz.




LB-363
Conf. 1 pz.



LB-364
Conf. 10 pz.



LB-365
Conf. 10 pz.



LB-366
Conf. 1 pz.



LB-367
Conf. 10 pz.



LB-368
Conf. 10 pz.



LB-369
Conf. 1 pz.


Scalpel Handles

Skalpellgriffe

Manici porta lame

Mangos de bisturí

Manches de bistouri



Für Mikro-Klingen
For Micro Blades
Per micro lame
Para micro Hojas
Pour Micro-Lame de bistouris




NK-1
100 mm




NK-3
130 mm



NK-2
150 mm



MA-B16
120 mm




GB-4
165 mm

Scalpel Blades and Handle

Skalpelle und Skalpellgriffe
Lame e manici per bisturi
Hojas y mangos de bisturí
Lames et manches de bistouri

Scalpels

Skalpelee
Bisturies
Bistouris
Bisturi



Packungen à 100
Packages of 100 ea
Scatole da 100 pezzi
Cajas de 100 c/u
Boîtes de 100 pièces




LB-10
Fig. 10



LB-11
Fig. 11



LB-12
Fig. 12



LB-12D
Fig. 12D



LB-13
Fig. 13



LB-15
Fig. 15



LB-15A
Fig. 15A



LB-15C
Fig. 15C



MA-B3



MA-B3M
Fig. 3M



MA-B7



MA-B7K




MA-B1
Fig. 5



MA-B2
Fig. 5A



MA-B8/MA-B9
155 mm



Doppelt Skalpellgriff
Scalpel Handle double Blade
Manico per bisturi a doppia lama
Mango de bisturí doble hojas
Manche pour bistouri double-lames


Scalpel Blades

Skalpelle

Lame per bisturi

Hojas de bisturí

Lames de bistouri



No. 4

Packungen à 100
Packages of 100 ea
Scatole da 100 pezzi
Cajas de 100 c/u
Boîtes de 100 pièces



LB-18
Fig. 18




LB-20
Fig. 20



LB-21
Fig. 21



LB-23
Fig. 23



**PAQUETTE
BLADE HOLDER**
GB-2
160 mm


Palatal, Gingival,
Connective Tissue Grafts
Uses standard
double-edged
razor blade

Mucotomo




MA-B4
Fig. 4

Klingenentfemer
Scalpel blade removing forceps
Pinze togliilama
Pinza para quitar hojas de bisturíes
Pince pour enlever les lames de bistouris



GS-18
135 mm




GS-16
115 mm

Knives

Messer
Coltelli
Cuchillos
Couteaux


Scalpels

Skalpelee
Bisturies
Bistouris
Bisturi



Knives

Messer
Coltelli
Cuchillos
Couteaux



NK-30/NK-40



1/1

NK-30
120 mm



1/1

NK-32
125 mm



1/1

NK-34
130 mm



1/1

NK-36
120 mm



1/1

NK-38
125 mm



1/1


NK-40
130 mm

SCISSORS




Micro Scissors

Mikro-Scheren
Micro Forbici
Micro Tijeras
Micro Ciseaux




FD-10
120 mm
Spitz
Sharp
Acuta
Aguda
Pointu

NOYES
FD-10
120 mm




FD-12
120 mm
Spitz
Sharp
Acuta
Aguda
Pointu

NOYES
FD-12
120 mm




FD-41 **FD-42**
120 mm 120 mm
Stumpf
Blunt
Smussa
Romo
Mousse

FD-41/FD-42
120 mm




FD-15 **FD-13**
120 mm 120 mm
Spitz
Sharp
Acuta
Aguda
Pointu

WESTCOTT
FD-13/FD-15
120 mm



FD-64 **FD-65**
120 mm 120 mm
Stumpf
Blunt
Smussa
Romo
Mousse

WESTCOTT
FD-64/FD-65
120 mm



FD-137 **FD-138**
150 mm 150 mm


FD-137/FD-138
150 mm

Scissors


Scheren
Forbici
Tijeras
Ciseaux

Scissors


Scheren
Forbici
Tijeras
Ciseaux




FD-182/FD-183
160 mm




FD-182
160 mm
Spitz
Sharp
Acuta
Aguda
Pointu




FD-183
160 mm
Spitz
Sharp
Acuta
Aguda
Pointu



1/1



FD-224/FD-225
145 mm




FD-224
145 mm
Spitz
Sharp
Acuta
Aguda
Pointu



FD-225
145 mm
Spitz
Sharp
Acuta
Aguda
Pointu


1/1




FD-145
160 mm

Gezahnt
Serrated
Dentellata
Dentada
Dentelé

Spitz
Sharp
Acuta
Aguda
Pointu




FD-06
160 mm




1/1

FD-06
160 mm

Spitz
Sharp
Acuta
Aguda
Pointu



FD-08
160 mm





1/1

FD-08
160 mm

Spitz
Sharp
Acuta
Aguda
Pointu

Scissors

Scheren
Forbici
Tijeras
Ciseaux





Scissors

Scheren
Forbici
Tijeras
Ciseaux

Scissors

Scheren
Forbici
Tijeras
Ciseaux


L = links
left
mancini
izquierdo
gauche



Scissors


Scheren
Forbici
Tijeras
Ciseaux

Gezahnt - Serrated
Dentellata - Dentelé - Dentadas




FG-156/FG-157
110 mm


Gezahnt - Serrated
Dentellata - Dentelé - Dentadas




FG-158/FG-159
110 mm




KILNER
FG-58
115 mm



JAMESON-WERBER
FG-123
130 mm



JAMESON
FG-124
145 mm




FG-3
115 mm

Scissors


Scheren
Forbici
Tijeras
Ciseaux

Scissors


Scheren
Forbici
Tijeras
Ciseaux




MICRO-IRIS
FG-13/FG-15
85 mm



FG-23/FG-24
110 mm




FG-104
100 mm




FG-103/FG-105
110 mm

Scissors


Scheren
Forbici
Tijeras
Ciseaux




METZENBAUM
FM-126/FM-127
145 mm




METZENBAUM
FM-12/FM-109
150 mm



METZENBAUM
FG-28
145 mm



METZENBAUM
FG-14
145 mm



JOSEPH
FG-147/FG-148
150 mm


Scissors

Scheren
Forbici
Tijeras
Ciseaux

Scissors


Scheren
Forbici
Tijeras
Ciseaux

Stumpf
Blunt
Smussa
Romas
Mousse




FG-101
110 mm

Stumpf
Blunt
Smussa
Romas
Mousse




FG-102
110 mm

Stumpf
Blunt
Smussa
Romas
Mousse




FG-21/FG-22
115 mm


Spitz
Sharp
Acute
Agudas
Pointu



FG-36/FG-37
115 mm





LOCKLIN
LK-1
115 mm



LOCKLIN MOD.
LK-3
110 mm

Scissors

Scheren
Forbici
Tijeras
Ciseaux




Scissors

Scheren
Forbici
Tijeras
Ciseaux


Scissors

Scheren
Forbici
Tijeras
Ciseaux

Gezahnt - Serrated
Dentellata - Dentadas
Dentelé




LOCKLIN
LK-2
160 mm




LOCKLIN
LK-4
160 mm

Gezahnt - Serrated
Dentellata - Dentadas
Dentelé




KELLY
FK-1/FK-2
160 mm




KELLY
FK-87/FK-88
180 mm

Scissors

Scheren
Forbici
Tijeras
Ciseaux




STEVENS
FG-50/FG-53
110 mm



FG-50

Stumpf
Blunt
Smussa
Romas
Mousse


Gerade
Straight
Diritta
Rectas
Droites



FG-51

Stumpf
Blunt
Smussa
Romas
Mousse


Gebogen
Curved
Curva
Curvas
Courbes



FG-52

Spitz
Sharp
Acuta
Agudas
Pointu


Gerade
Straight
Diritta
Rectas
Droites




FG-53

Spitz
Sharp
Acuta
Agudas
Pointu

Gebogen
Curved
Curva
Curvas
Courbes




REYNOLDS
FG-54/FG-55



FG-54
150 mm


Gerade
Straight
Diritta
Rectas
Droites




FG-55
150 mm

Gebogen
Curved
Curva
Curvas
Courbes

Stumpf
Blunt
Smussa
Romas
Mousse




KILNER
FG-19C/FG-67




FG-66
120 mm

Gerade
Straight
Diritta
Rectas
Droites




FG-67
120 mm

Gebogen
Curved
Curva
Curvas
Courbes



FG-19R
135 mm

Gerade
Straight
Diritta
Rectas
Droites



FG-19C
135 mm

Gebogen
Curved
Curva
Curvas
Courbes


Scissors

Scheren
Forbici
Tijeras
Ciseaux


Scissors

Scheren
Forbici
Tijeras
Ciseaux

Gezahnt - Serrated
Dentellata - Dentadas
Dentelé




LA GRANGE
LG-1
115 mm




LA GRANGE
LG-6
155 mm


Gezahnt - Serrated
Dentellata - Dentada
Dentée



GOLDMAN-FOX S FORM
GFS
130 mm




NEUMANN
FG-20
130 mm




DEAN
DEAN
180 mm

Scissors


Scheren
Forbici
Tijeras
Ciseaux




QUINBY
FG-38
130 mm




CHADWICK
FG-39
115 mm




MAYO-LEXER
FM-79
165 mm



LEXER
FM-180/FM-181
165 mm



FM-180



FM-181

Scissors

Scheren
Forbici
Tijeras
Ciseaux


Scissors

Scheren
Forbici
Tijeras
Ciseaux



FM-9/FM-10
145 mm

FM-27/FM-63



MAYO
FM-5/FM-6
140 mm

FM-7/FM-42

Scissors

Scheren
Forbici
Tijeras
Ciseaux

SUPER CUT



Super Cut Scissors have specially designed cutting edges.
One regular edge and **one sharp-knife edge**.




Extremely Sharp Razor Edge / Atraumatic
Extrem Scharfe Spitzer Beck / Atraumatische
Tranchant extrêmement aiguisée / Atraumatique
Borde extremadamente afilado / Atraumática
Taglio estremamente affilato / Atraumatica





One side serrated
Eine seite gezahnt
Un côté dentelé
Un lado dentado
Un lato dentellato

Super Cut





METZENBAUM
FM-318/FM-319
180 mm

Super Cut



METZENBAUM
FM-298/FM-299
125 mm

Super Cut




STEVENS
FM-300/FM-301
125 mm

Scissors


Scheren
Forbici
Tijeras
Ciseaux

Scissors

Scheren
Forbici
Tijeras
Ciseaux




BABY - METZENBAUM
FM-47/FM-48
115 mm




TOENNIS-ADSON
FG-60/FG-61
180 mm

Präparierscheren - Dissecting Scissors
Forbici per dissezione - Tijeras de disección - Ciseaux à dissection





METZENBAUM
FM-1/FM-2
145 mm



METZENBAUM-FINO
FM-1F/FM-2F
145 mm

Scissors

Scheren
Forbici
Tijeras
Ciseaux




Scissors


Scheren
Forbici
Tijeras
Ciseaux

Scissors


Scheren
Forbici
Tijeras
Ciseaux




METZENBAUM-FINO
FM-201/FM-204




FM-201
140 mm



FM-202
140 mm




FM-203
180 mm




FM-204
180 mm




Tungsten-Carbide Inserts
Mit Hartmetalleinsätzen
Plaquettes en carbure de tungstène
Plaquetas de carburo de tungsteno
Inserti in carburo di tungsteno




METZENBAUM-FINO
FM-96/FM-99




FM-96
140 mm



FM-97
180 mm



FM-98
140 mm



FM-99
180 mm

Scissors


Scheren
Forbici
Tijeras
Ciseaux




Scheren
Scissors
Forbici
Tijeras
Ciseaux




Tungsten-Carbide Inserts
Mit Hartmetalleinsätzen
Plaquettes en carbure de tungstène
Plaquetas de carburo de tungsteno
Inserti in carburo di tungsteno




BABY-METZENBAUM
FM-118
115 mm




BABY-METZENBAUM
FM-119
115 mm



KILNER
FM-350
125 mm




IRIS
FG-1TC/FG-2TC
115 mm



FG-167/FG-168
110 mm

Gebogen - Curved
Curva - Curvas
Courbes



GREENBERG-PAR
FM-317
115 mm


Gezahnt
Serrated
Dentellata
Dentadas
Dentelé

Scissors

Scheren
Forbici
Tijeras
Ciseaux


Scissors

Scheren
Forbici
Tijeras
Ciseaux



LK-10
160 mm

One side serrated
Eine Seite gezahnt
Un côté dentelé
Un lado dentado
Un lato dentellato




LK-12
160 mm

Non Serrated
Nicht Gezahnt
Non dentelé
No dentado
No dentellato




Eine Schneide gezahnt S-förmig gebogen
One toothed cutting edge S-shaped
Una lama dentata curva a forma S
Una hoja dentada forma S
Une lame dentelée courbée en S




**KELLY
FK-1SC/FK-2SC**
160 mm

Eine Schneide gezahnt
One toothed cutting
Una lama dentata
Una hoja dentada
Une lame dentelé



**KELLY
FK-1TC/FK-2TC**
160 mm

Eine Schneide gezahnt
One toothed cutting
Una lama dentata
Una hoja dentada
Une lame dentelé




DEAN-TC
175 mm

Scissors

Scheren
Forbici
Tijeras
Ciseaux



Wellenschliff
TC




IRIS
FM-100L
110 mm

Links Scheren
Left Scissors
Forbici per mancini
Tijeras para zurdos
Gauche Ciseaux




Tungsten-Carbide Inserts
Mit Hartmetalleinsätzen
Plaquettes en carbure de tungstène
Plaquetas de carburo de tungsteno
Inserti in carburo di tungsteno

Wellenschliff
TC




METZENBAUM
FM-102
145 mm

Wellenschliff
TC



METZENBAUM
FM-104
185 mm

Wellenschliff
TC



IRIS
FM-100
110 mm

Scissors


Scheren
Forbici
Tijeras
Ciseaux

Scissors


Scheren
Forbici
Tijeras
Ciseaux



Extremely Sharp Razor Edge / Atraumatic
Extrem Scharfe Spitzer Beck / Atraumatische
Tranchant extrêmement aiguës / Atraumatique
Borde extremadamente afilado / Atraumática
Taglio estremamente affilato / Atraumatica



FG-1SCZ/FG-2SCZ
110 mm



IRIS
FG-1SC/FG-2SC
110 mm

METZENBAUM
FM-1SC/FM-2SC
145 mm

MAYO
FM-5SC/FM-6SC
140 mm

METZENBAUM
FM-3SC/FM-4SC
185 mm

Scissors

Scheren
Forbici
Tijeras
Ciseaux



SUPER CUT



Super Cut Scissors have specially designed cutting edges.
One regular edge and **one sharp-knife edge**.






Extremely Sharp Razor Edge / Atraumatic

Extrem Scharfe Spitzer Beck / Atraumatische
Tranchant extrêmement aiguisée / Atraumatique
Borde extremadamente afilado / Atraumática
Taglio estremamente affilato / Atraumatica



One side serrated

Eine seite gezahnt
Un côté dentelé
Un lado dentado
Un lato dentellato




Scissors



Scheren
Forbici
Tijeras
Ciseaux

Scissors

Scheren
Forbici
Tijeras
Ciseaux




JOSEPH
FG-153/FG-154
150 mm




GOLDMAN-FOX
GFC-TC/GFR-TC
135 mm

GFR-TC GFC-TC
Gezahnt - Serrated
Dentellata - Dentadas
Dentelé



PECK-JOSEPH
FG-92/FG-98
140 mm


FG-92 FG-98
Spitz
Sharp
Acuta
Aguda
Pointu
Stumpf
Blunt
Smussa
Romas
Mousse




GREGORY
FM-339
140 mm

Scissors

Scheren
Forbici
Tijeras
Ciseaux



FM-146
145 mm




Ligature Scissors


Ligatur Scheren
Forbici per legature
Tijeras para ligaduras
Ciseaux à ligatures

Scissors


Scheren
Forbici
Tijeras
Ciseaux




SPENCER
SC-9
90 mm




SPENCER
SC-1
110 mm




SPENCER
SC-5
110 mm




SPENCER
SC-2
145 mm




SPENCER
SC-14
120 mm



LITTAUER
SC-15
140 mm



SPENCER
SC-13
130 mm



KELLY
SC-16
150 mm


Wire and plate cutting scissors

Draht und blechscheren

Forbici per fili metallici


Tijeras para cortar alambre

Ciseaux coupe-fils métalliques



Draht-und Blechschere gezahnt
Wire cutting scissor serrated
Forbice universale dentellata
Tijeras para cortar hilo de acero dentadas
Coupe-fils dentelé

DP-90
120 mm



GS-26
110 mm

Spitz
Sharp
Acute
Agudas
Pointu




BEEBE
DP-91/DP-101
110 mm

Stumpf
Blunt
Smusse
Romas
Mousse




BEEBE
DP-102/DP-105
110 mm

TWEEZERS




Delicate Forceps


Mikro Pinzetten
Micro Pinzette
Pinzas finas
Micro Pincers




BD-360
110 mm




BD-361
110 mm




BD-362
100 mm




BD-363
100 mm




BD-364
110 mm



BD-365
110 mm



BD-366
130 mm



BD-367
130 mm


Delicate Forceps

Mikro Pinzetten
Micro Pinzette
Pinzas finas
Micro Pincers

Tweezers


Pinzette
Pince À Épiler
Pinzette
Pintset

Tissue Forceps
1x2 Teeth



BISHOP-HARMON
MF-177 1,2 mm
85 mm

Dressing Forceps
criss-cross serrated tips



BISHOP-HARMON
MF-249 1,2 mm
85 mm



Diamond Tip
Diamantiert
Punta diamantata
Punta de diamante
Pointe de diamant

1/1



1x2 Zähne
1x2 Teeth
1x2 Denti
1x2 Dientes
1x2 Dents

CASTROVIEJO
BD-348 0.9 mm
120 mm


1/1



9x9 Zähne
9x9 Teeth
9x9 Denti
9x9 Dientes
9x9 Dents


BROWN
BD-349
120 mm

1,2 mm




ADSON
BD-207
120 mm

1,2 mm




ADSON
BD-208
120 mm

1/1



ADSON-BROWN
BD-209
120 mm

1/1





ADSON-BROWN
BD-253
120 mm

TC




Delicate Forceps

Mikro Pinzetten
Micro Pinzette
Pinzas finas
Micro Pincers




MF-110/MF-111
120 mm




MF-110

1/1




MF-111

1/1




MF-112/MF-113
150 mm




MF-112

1/1




MF-113

1/1




MF-114/MF-115
180 mm



MF-114

1/1



MF-115

1/1

Micro Forceps

Mikro Pinzetten
Micro Pinzette
Micro Pinzas
Micro Pincers

Tweezers


Pinzette
Pince À Épiler
Pinzette
Pintset

Mikro Pinzette - Rundgriff Ø 8 mm
Mit Führungsstift diamantiert
Fadenplatte 6x0,7 mm

Micro Forceps
Micro Pinzette
Micro Pinzas
Micro Pincers



Diamond Tip
Diamantiert
Punta diamantata
Punta de diamante
Pointe de diamant



MF-370/MF-371
150 mm

MF-372/MF-373
180 mm

MF-374/MF-375
210 mm

MF-376/MF-377
230 mm

Micro Forceps

Mikro Pinzetten
Micro Pinzette
Micro Pinzas
Micro Pince



Mikro Pinzette - Rundgriff Ø 8 mm
Mit Führungsstift diamantiert
Fadenplatte 6x0,7 mm



Diamond Tip
Diamantiert
Punta diamantata
Punta de diamante
Pointe de diamant

1x2 Zähne
1x2 Teeths
1x2 Denti
1x2 Dientes
1x2 Dents



MF-378/MF-379
150 mm

1x2 Zähne
1x2 Teeths
1x2 Denti
1x2 Dientes
1x2 Dents



MF-380/MF-381
180 mm

1x2 Zähne
1x2 Teeths
1x2 Denti
1x2 Dientes
1x2 Dents



MF-382/MF-383
210 mm


1x2 Zähne
1x2 Teeths
1x2 Denti
1x2 Dientes
1x2 Dents




MF-384/MF-385
230 mm




MF-384 MF-385



MF-382 MF-383



MF-380 MF-381




MF-378 MF-379

Titanium Tissue Forceps


Titan-Pinzetten
Pinzette in titanio
Pinzas de titanio
Pinces en titane

Tweezers


Pinzette
Pince À Épiler
Pinzette
Pintset




BD-285T
90 mm




BD-270T
130 mm




ADSON
BD-30T
120 mm




ADSON
BD-31T
120 mm




BD-09T
130 mm



GERALD
MF-101T/MF-102T
180 mm



DE BAKEY
BD-80T/BD-160T



COLLEGE
TT-01
150 mm


BD-75T
145 mm

Delicate Forceps

Mikro Pinzetten
Micro Pinzette
Pinzas finas
Micro Pincers




Rund-Griffe Ø 8 mm
Round handle Ø 8 mm
Manico rotondo Ø 8 mm
Mango redondo Ø 8 mm
Manche rond Ø 8 mm




MF-124 - 150 mm
MF-125 - 180 mm
MF-126 - 210 mm

MF-124/MF-129



MF-127 - 150 mm
MF-128 - 180 mm
MF-129 - 210 mm


Titan
Titanium
Titanio
Titanio
Titane



MF-142T - 150 mm
MF-143T - 180 mm
MF-144T - 210 mm

Ø 2 mm

MF-145T - 150 mm
MF-146T - 180 mm
MF-147T - 210 mm




MF-142T/MF-147T

Cotton and Dressing Pliers

Zahnpinzetten
Pinzette per medicazione
Pinzas para vendajes
Précettes pour cotton

Tweezers

Pinzette
Pince À Épiler
Pinzette
Pintset




Endodontic Forceps


Pinzette für Endodontie
Pinzette per endodonzia
Pinzas para endodoncia
Pinces pour endodontie




Endo Forceps



1/1



1/1




1/1




1/1



DA-04
Endo
150 mm



DA-05
Endo / Locking
150 mm



DA-09
Endo / Locking
150 mm




DA-11
Endo / Locking
155 mm

Tissue and Dressing Forceps

Chirurgische und Anatomische Pinzetten
Pinzette anatomiche e chirurgiche
Pinzas de diseccion y quirùrgicas
Pincas anatomiques chirurgicales


Tweezers

Pinzette
Pince À Épiler
Pinzette
Pintset



Atraumatische pinzetten
Non-traumatic forceps
Pinzette atraumatiche
Pinzas atraumaticas
Pincas atraumatiques

DE BAKEY
BD-80/BD-165



Tissue Forceps

Chirurgische Pinzetten
Pinzette chirurgiche
Pinzas quirùrgicas
Pincas quirurgicales




Taschenmarkierungspinzetten

Pocket markers

Pinzette segnatasche

Pinzas para marcar el sulco periodontal


Pincas à marquer la poche



BD-279
150 mm

BD-283
150 mm


McINDOE
BD-279/BD-283
150 mm



BD-280
150 mm

BD-284
150 mm

GILLIES
BD-280/BD-284
150 mm



CRANE-CAPLAN
DB-01
150 mm






CRANE-CAPLAN
DB-02
150 mm

Tissue Forceps

Chirurgische Pinzetten
Pinzette chirurgiche
Pinzas quirùrgicas
Pincas quirurgicales

Tweezers

Pinzette
Pince À Épiler
Pinzette
Pintset



Dressing Forceps

Anatomische Pinzetten


Pinzette anatomiche

Pinzas de disección

Pinces anatomiques




BD-270
130 mm



1/1




BD-145
130 mm




1/1



BD-16/BD-148
160 mm



1/1




1/1

BD-16
160 mm

BD-148
160 mm



BD-10/BD-76



Mod. Fine


BD-10
145 mm

BD-75
145 mm

BD-76
200 mm



BD-97/BD-99



Standard

BD-97
100 mm

BD-98
115 mm


BD-99
130 mm

Tissue Forceps


Chirurgische Pinzetten
Pinzette chirurgiche
Pinzas quirúrgicas
Pincas quirúrgicas

Tweezers


Pinzette
Pince À Épiler
Pinzette
Pintset




BD-240/BD-249



BD-60/BD-62
110 mm



STILLE
BD-61/BD-67



BD-142/BD-143


Tissue Forceps

Chirurgische Pinzetten


Pinzette chirurgiche

Pinzas quirúrgicas


Pinces cirurgicales




BD-40
110 mm




BD-41
110 mm



BD-42/BD-43
110 mm




SEMKEN
BD-44/BD-45
130 mm



SEMKEN
BD-46/BD-47
130 mm

1x2 Zähne
1x2 Teeth
1x2 Denti
1x2 Dents
1x2 Dientes



SEMKEN
BD-91/BD-95
150 mm


1x2 Zähne
1x2 Teeth
1x2 Denti
1x2 Dents
1x2 Dientes

Tissue Forceps


Chirurgische Pinzetten
Pinzette chirurgiche
Pinzas quirúrgicas
Pincas quirúrgicas

Tweezers


Pinzette
Pince À Épiler
Pinzette
Pintset




BD-109/BD-144



BD-21/BD-282




BD-38/BD-134
200 mm




BD-19/BD-137

Tissue Forceps

Chirurgische Pinzetten
Pinzette chirurgiche
Pinzas quirúrgicas
Pincers chirurgicales




GERALD
MF-101
180 mm




GERALD
MF-102
180 mm



GERALD
MF-103
180 mm



GERALD
MF-104
180 mm







GERALD-BROWN
MF-105
180 mm

Tissue Forceps


Chirurgische Pinzetten
Pinzette chirurgiche
Pinzas quirúrgicas
Pincers chirurgicales

Tweezers

Pinzette
Pince À Épiler
Pinzette
Pintset




Pinzetten zum Halten Membranen
Membrane Holding forceps
Pinza per fissaggio membrane
Pinza para sostener membranas
Pince pour le maintien membranes






Tissue Forceps


Chirurgische Pinzetten
Pinzette chirurgiche
Pinzas quirúrgicas
Pincas quirúrgicas




ADSON-BROWN
BD-92/BD-141
120 mm




ADSON
BD-90
120 mm




BROWN
BD-94
150 mm



ADSON
BD-380
150 mm



ADSON
BD-381
150 mm




ADSON-BROWN
BD-382
150 mm

Tissue Forceps

Chirurgische Pinzetten
Pinzette chirurgiche
Pinzas quirúrgicas
Pincas chirurgicales


Tweezers

Pinzette
Pince À Épiler
Pinzette
Pintset




Mit Fadenplatte
With tying platform
Con piano prendifilo
Con plataforma
Avec plateaux pour monofil

ADSON
BD-370
120 mm




Mit Fadenplatte
With tying platform
Con piano prendifilo
Con plataforma
Avec plateaux pour monofil
1x2 Zähne
1x2 Teeth
1x2 Denti
1x2 Dientes
1x2 Dents

ADSON
BD-371
120 mm



ADSON
BD-30/BD-73



ADSON
BD-31/BD-72


Tissue Forceps - Dressing forceps

Chirurgische Pinzetten - Anatomische Pinzetten

Pinzette chirurgiche - Pinzette anatomiche

Pinzas quirúrgicas - Pinzas de disección


Pinces chirurgicales - Pinces anatomiques



Extra fein
Extra fine
Extra fine
Extra finas
Extra fins

Gebogen
Curved
Curva
Curvas
Courbé

**MICRO-ADSON
BD-372**
120 mm




Extra fein
Extra fine
Extra fine
Extra finas
Extra fins

Gebogen
Curved
Curva
Curvas
Courbé


1x2 Zähne
1x2 Teeth
1x2 Denti
1x2 Dientes
1x2 Dents

**MICRO-ADSON
BD-373**
120 mm




2x3 Zähne
2x3 Teeth
2x3 Denti
2x3 Dientes
2x3 Dents


**ADSON
BD-377**
120 mm



**BEASLEY-BABCOCK
BD-51**
120 mm




**GRAEFE
BD-135**
110 mm



Gerieft
Serrated
Zigrinata
Estriada
Striée

Abgewinkelt
Angled
Angolata
Angulada
Coudée en haut

**ADSON
BD-378**
120 mm



1x2 Zähne
1x2 Teeth
1x2 Denti
1x2 Dientes
1x2 Dents

Abgewinkelt
Angled
Angolata
Angulada
Coudée en haut

**ADSON
BD-379**
120 mm

Tissue Forceps - Dressing forceps

Chirurgische Pinzetten - Anatomische Pinzetten
Pinzette chirurgiche - Pinzette anatomiche
Pinzas quirùrgicas - Pinzas de diseccion
Pincas chirurgicales - Pincas anatomiques

Tweezers

Pinzette
Pince À Épiler
Pinzette
Pintset



BD-09
130 mm

BD-18
180 mm

BD-20
200 mm

BD-09/BD-20



Pinzas quirùrgicas - Pinzas de diseccion
Pincas chirurgicales - Pincas anatomiques

BD-23
230 mm

BD-25
250 mm


BD-26
300 mm

BD-23/BD-26



1x2 Zähne
1x2 Teeth
1x2 Denti
1x2 Dientes
1x2 Dents

BD-271/BD-276




Tissue Forceps

Chirurgische Pinzetten


Pinzette chirurgiche

Pinzas quirúrgicas


Pinces cirurgicales




BD-290/BD-292



BD-351/BD-356



**POTTS-SMITH
BD-14/BD-53**
180 mm










**POTTS-SMITH
BD-287/BD-289**

Tissue Forceps, Non-traumatic Forceps

Chirurgische Pinzetten, Atraumatische Pinzetten
Pinzette chirurgische, Pinzette atraumatice
Pinzas quirùrgicas, Pinzas atraumàticas
Pincas chirurgicales, Pincas atraumatiques

Tweezers

Pinzette
Pince À Épiler
Pinzette
Pintset



Tissue Forceps, Non-traumatic Forceps

Chirurgische Pinzetten, Atraumatische Pinzetten

Pinzette chirurgiche, Pinzette atraumatiche


Pinzas quirùrgicas, Pinzas atraumàticas

Pinces chirurgicales, Pinces atraumatiques




BD-210/BD-213

Standard




SEMKEN
BD-216/BD-217



Tungsten-Carbide Inserts
Mit Hartmetalleinsätzen
Plaquettes en carbure de tungstène
Plaquetas de carburo de tungsteno
Inserti in carburo di tungsteno



BD-221/BD-224



DE BAKEY
BD-267
160 mm



Tissue Forceps, Non-traumatic Forceps

Chirurgische Pinzetten, Atraumatische Pinzetten

Pinzette chirurgiche, Pinzette atraumatiche

Pinzas quirùrgicas, Pinzas atraumàticas

Pinces chirurgicales, Pinces atraumatiques

Tweezers


Pinzette

Pince À Épiler


Pinzette

Pintset

Naht-und Membranpinzette - Suture and membrane Forceps Pinzette per suture e membrane - Pinzas para suturas y membranas Pincettes à suture et à membrane



Artikulationspapier-Pinzetten - Articulating Paper Forceps Pinza per carta articolazione - Pinzas para papel de articular Pince pour papier à articuler



DRIVERS



Needle Holders


Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles

Drivers


Treiber
Conductrices
Autiste
Drayvery

Mit Scheren
With Scissor
Con forbice
Con tijeras
Avec ciseaux


Glattes Maul
Smooth Jaw
Morso liscio
Mandibula lisa
Mors lisse




STEVENS
FD-38
160 mm



FD-26
140 mm



FD-14/FD-16
140 mm



CASTROVIEJO
FD-17/FD-19
140 mm

Needle Holders


Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles




Mit Locking-System
With Locking System (stop)
Con sistema di bloccaggio
Con sistema de bloqueo
Avec système de verrouillage

Glattes Maul
Smooth Jaw
Morso liscio
Mandíbula lisa
Mors lisse


Feingerieftes Maul
Fine serrated jaws
Morso zigrinato fine
Mandíbula estriada fina
Mors strié fin




**BARRAQUER
FD-216**
150 mm



**BARRAQUER
FD-217**
150 mm



**BARRAQUER
FD-236/FD-237**
150 mm



**BARRAQUER
FD-24**
150 mm

Needle Holders


Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles

Drivers

Treiber
Conductrices
Autiste
Drayvery



Tungsten-Carbide Inserts
Mit Hartmetalleinsätzen
Plaquettes en carbure de tungstène
Plaquetas de carburo de tungsteno
Inserti in carburo di tungsteno



Needle Holders

Nadelhalter

Port'aghi




Portaagujas

Porte aiguilles



Tungsten-Carbide Inserts

Mit Hartmetalleinsätzen
Plaquettes en carbure de tungstène
Plaquetas de carburo de tungsteno
Inserti in carburo di tungsteno



Needle Holders




Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles

Drivers

Treiber
Conductrices
Autiste
Drayvery



Tungsten-Carbide Inserts
Mit Hartmetalleinsätzen
Plaquettes en carbure de tungstène
Plaquetas de carburo de tungsteno
Inserti in carburo di tungsteno




Needle Holders

Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles




Tungsten-Carbide Inserts
Mit Hartmetalleinsätzen
Plaquettes en carbure de tungstène
Plaquetas de carburo de tungsteno
Inserti in carburo di tungsteno


Glattes Maul
Smooth jaws
Con morso liscio
Mandibula lisa
A' mors lisse



CASTROVIEJO
FD-21/FD-23
140 mm



MICRO-RYDER
PH-140/PH-142



DIETRICH
FD-33/FD-35

Micro Needle Holders - Micro Scissors

Mikro Nadelhalter - Mikro Scheren
 Micro Port'aghi - Micro Forbici
 Micro Portaagujas - Micro Tijeras
 Micro Porte aiguilles - Micro Ciseaux

Drivers


Treiber
 Conductrices
 Autiste
 Drayvery



Rund-Griffe Ø 8 mm
 Round handle Ø 8 mm
 Manico rotondo Ø 8 mm
 Mango redondo Ø 8 mm
 Manche rond Ø 8 mm


Mikro Nadelhalter
 Micro Needle Holders
 Micro port'aghi
 Micro portaagujas
 Micro porte aiguilles

Mikro Scheren
 Micro Scissors
 Micro forbici
 Micro tijeras
 Micro ciseaux



FD-562T FD-563T


FD-562T/FD-563T
 180 mm



FD-574T FD-575T

FD-574T/FD-575T
 180 mm

Mit Sperre
 With Lock
 Con arresto
 Con cierre
 Avec Arrêt



FD-587T FD-588T

FD-587T/FD-588T
 180 mm

Needle Holders

Nadelhalter

Port'aghi

Portaagujas

Porte aiguilles



Tungsten-Carbide Inserts

Mit Hartmetalleinsätzen

Plaquettes en carbure de tungstène

Plaquetas de carburo de tungsteno

Inserti in carburo di tungsteno




Diamond Tip

Diamantiert

Punta diamantata

Punta de diamante

Pointe de diamant




FD-730
150 mm

FD-731
170 mm


CASTROVIEJO
Fine
Straight
Locking
Thick Jaw
Round Handle

CASTROVIEJO
FD-730 / FD-731
150 / 170 mm



FD-599 FD-600

FD-599/FD-600
150 mm



FD-601 FD-602


FD-601/FD-602
180 mm

Needle Holders


Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles

Drivers

Treiber
Conductrices
Autiste
Drayvery




DERF
PH-39
125 mm




PH-38/PH-64
120 mm

Mit Glatten maul
With smooth jaws
Con morso liscio
Con mandibula lisa
Avec mors lisse




HALSEY
PH-48/PH-63
130 mm




HALSEY
PH-37
130 mm

Mit Glatten maul
With smooth jaws
Con morso liscio
Con mandibula lisa
Avec mors lisse




CRILE-WOOD
PH-61
150 mm




CRILE-MURRAY
PH-62
150 mm

Needle Holders

Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles




MINI-OLSEN-HEGAR
BM-61
120 mm




OLSEN-HEGAR
BM-62
140 mm

Nadelhalter mit Schere - Needle holders with scissors
Port'aghi con forbici - Portaagujas con tijeras
Porte aiguilles avec ciseaux



OLSEN-HEGAR
BM-63
170 mm



MINI-GILLIES
PO-59
145 mm


Nadelhalter mit Schere - Needle holders with scissors
Port'aghi con forbici - Portaagujas con tijeras
Porte aiguilles avec ciseaux

Needle Holders

Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles

Drivers


Treiber
Conductrices
Autiste
Drayvery




RYDER
PH-80/PH-82



Tungsten-Carbide Inserts
Mit Hartmetalleinsätzen
Plaquettes en carbure de tungstène
Plaquetas de carburo de tungsteno
Inserti in carburo di tungsteno



RYDER
PH-83
140 mm




PH-84
140 mm

Needle Holders


Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles




Tungsten-Carbide Inserts
Mit Hartmetalleinsätzen
Plaquettes en carbure de tungstène
Plaquetas de carburo de tungsteno
Inserti in carburo di tungsteno



CRILE-WOOD
PH-79
145 mm



FINE SWEDISH PATTERN
PH-85
150 mm



SWEDISH MOD.
PH-106
150 mm

Needle Holders


Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles

Drivers


Treiber
Conductrices
Autiste
Drayvery




Tungsten-Carbide Inserts
Mit Hartmetalleinsätzen
Plaquettes en carbure de tungstène
Plaquetas de carburo de tungsteno
Inserti in carburo di tungsteno




PH-53/PH-65
120 mm




**WEBSTER
PH-66**
130 mm




**HALSEY
PH-67**
130 mm



**DERF
PH-69**
115 mm



**CONVERSE
PH-68**
130 mm




**HALSEY
PH-73**
130 mm

Needle Holders


Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles



Tungsten-Carbide Inserts
Mit Hartmetalleinsätzen
Plaquettes en carbure de tungstène
Plaquetas de carburo de tungsteno
Inserti in carburo di tungsteno




DE BAKEY
PH-110/PH-111




MICRO-VASCULAR
PH-101/PH-102/PH-103

Vascular
Vaskuläre
Vasculaire
Vasculares
Vascolari



HEGAR
PH-120/PH-121/PH-122



MAYO-HEGAR-FINO
PH-127/PH-128/PH-129

Needle Holders


Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles

Drivers


Treiber
Conductrices
Autiste
Drayvery




Tungsten-Carbide Inserts
Mit Hartmetalleinsätzen
Plaquettes en carbure de tungstène
Plaquetas de carburo de tungsteno
Inserti in carburo di tungsteno




HEGAR
PH-59/PH-72/PH-75/PH-76



PH-58
150 mm



SAROT
PH-132/PH-133/PH-134




CRILE-WOOD
PH-77
180 mm

Needle Holders


Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles




Tungsten-Carbide Inserts
Mit Hartmetalleinsätzen
Plaquettes en carbure de tungstène
Plaquetas de carburo de tungsteno
Inserti in carburo di tungsteno




PH-99
120 mm




**CRILE-WOOD
PH-100**
150 mm



**BABY-CRILE-WOOD
PH-54**
150 mm



**CRILE-WOOD
PH-71**
150 mm



**CRILE-WOOD
PH-90**
150 mm

Needle Holders


Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles

Drivers


Treiber
Conductrices
Autiste
Drayvery




Tungsten-Carbide Inserts
Mit Hartmetalleinsätzen
Plaquettes en carbure de tungstène
Plaquetas de carburo de tungsteno
Inserti in carburo di tungsteno




MINI-GILLIES
PO-99
145 mm




MINI-GILLIES
PO-99L
145 mm



FOSTER-GILLIES
PO-104/PO-105
120 mm




OLSEN-HEGAR
BM-01/BM-02




OLSEN-HEGAR
BM-03/BM-04

Needle Holders


Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles




ARRUGA
PO-85/PO-87
160 mm




BOYNTON
PO-55/PO-56
120 mm




BOYNTON
PO-88
120 mm




BOYNTON
PO-94
125 mm



MATHIEU
PO-50/PO-52



MATHIEU
PO-82/PO-84



LICHTENBERG
PO-86
175 mm


Needle Holders

Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles


Drivers

Treiber
Conductrices
Autiste
Drayvery

Needle Holders
Port'aghi
Portaagujas
Porte aiguilles




CRILE
PO-37
150 mm




TÖNNIS
PO-38
190 mm

Schlankesmodell - Delicate pattern - Modello fine
Modelo delegado - Modèle fin




MAYO-HEGAR
PH-70/PH-74




MAYO-HEGAR
PH-40/PH-42

Needle Holders


Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles




MATHIEU
PO-65/PO-81



MATHIEU
PO-79/PO-80
140 mm



MATHIEU
PO-46/PO-48




PO-40
140 mm

Needle Holders


Nadelhalter
Port'aghi
Portaagujas
Porte aiguilles

Drivers


Treiber
Conductrices
Autiste
Drayvery




MATHIEU
PO-45
125 mm




MATHIEU
PO-60/PO-71
140 mm



SMAHA
PO-68
140 mm



SMAHA
PO-69
140 mm



SMAHA
PO-76/PO-77
140 mm

GRASPING FORCEPS



Towel Clamps

Tuchklemmen

Fissateli

Pinzas de campo

Pinces fixe-compresses


Hemostats /Grasping Forceps

Hämostatika / Greifen


Hémostats / Saisir

Emostatici / Presa


Krovoostanavlivayushchiye




SCHAEFFEL
BF-38
85 mm




JONES
BF-48
90 mm




FAURE
BF-40
85 mm



BACKHAUS
BF-42/BF-46




BF-42
90 mm




BF-41
100 mm




BF-44
110 mm



BF-46
130 mm



BF-43
135 mm



BF-45
140 mm

Towel Clamps

Tuchklemmen


Fissateli

Pinzas de campo


Pinces fixe-compresses




Für Papiertücher
For paper cloths
Per teli di carta
Para paños
Pour papier



LORNA
BF-50
100 mm



BF-52
130 mm




BF-75
140 mm

Sponge and Dressing Forceps


Schwamm und Tamponzangen
Pinze per tamponi e sponge
Pinzas para esponja y tampones
Pincers pour éponges et à pansement

Hemostats /Grasping Forceps


Hämostatika / Greifen
Hémostats / Saisir
Emostatici / Presa
Krovoostanavlivayushchiye




PEAN-MURPHY
BF-22/BF-24



PEAN-MURPHY
BF-23/BF-25




GROSS
BF-94/BF-95
200 mm



GROSS
BF-10 130 mm
BF-12 160 mm

Sinus klemmen
Sinus forceps
Pinze Sinus
Pinzas para seno
Pincers à Sinus



LISTER
BF-20
150 mm


Hemostatic Forceps

Arterienklemmen


Pinze emostatiche

Pinzas hemostáticas


Pinces hémostatiques




HARTMANN/BABY-MOSQUITO
P-120/P-121
95 mm




HARTMANN/BABY-MOSQUITO
P-264/P-265
95 mm



DANDY-MOSQUITO
P-332
120 mm



JACOBSON-MICRO
P-266
130 mm




JACOBSON-MICRO
P-267
130 mm

Hemostatic Forceps


Arterienklemmen
Pinze emostatiche
Pinzas hemostáticas
Pincas hémostatiques

Hemostats /Grasping Forceps


Hämostatika / Greifen
Hémostats / Saisir
Emostatici / Presa
Krovoostanavlivayushchiye




P-15
130 mm



**PEET
DR-40**
120 mm



**MICRO-HALSTED
P-101/P-102**



**MICRO-HALSTED
P-135**
115 mm


Hemostatic Forceps

Arterienklemmen


Pinze emostatiche

Pinzas hemostáticas


Pinces hémostatiques




HALSTED
P-91/P-92
125 mm



HALSTED
P-93/P-94
125 mm



HALSTED-MOSQUITO
P-96/P-97
120 mm




HALSTED-MOSQUITO
P-95/P-98
120 mm

Hemostatic Forceps


Arterienklemmen
Pinze emostatiche
Pinzas hemostáticas
Pincas hémostatiques

Hemostats /Grasping Forceps

Hämostatika / Greifen
Hémostats / Saisir
Emostatici / Presa
Krovoostanavlivayushchiye




P-215/P-217
160 mm




P-215
160 mm

P-216
160 mm

P-217
160 mm




BABY-MIXTER
P-136/P-137




P-136
140 mm

P-137
180 mm



PEAN
P-250/P-253




P-250
130 mm


P-252
140 mm

P-251
130 mm

P-253
140 mm




PEAN
P-256/P-257




P-256
145 mm

P-257
145 mm



KOCHER
P-260/P-261



P-260
145 mm

P-261
145 mm


Hemostatic Forceps

Arterienklemmen


Pinze emostatiche

Pinzas hemostáticas


Pinces hémostatiques




CRILE
P-81/P-82
140 mm



CRILE
P-86/P-87
160 mm



KELLY
P-84/P-85
140 mm




KOCHER
P-49/P-89

Hemostatic Forceps


Arterienklemmen
Pinze emostatiche
Pinzas hemostáticas
Pincas hémostatiques

Hemostats /Grasping Forceps


Hämostatika / Greifen
Hémostats / Saisir
Emostatici / Presa
Krovoostanavlivayushchiye




LERICHE
P-61/P-62
150 mm



LERICHE
P-63/P-64
150 mm



BABY-CRILE
P-230/P-233
140 mm



HALSTED
P-238/P-241
185 mm



Tissue Grasping Forceps

Gewebefasszangen

Pinze per tessuti

Pinzas para tejidos

Pinces pour tissus




RETRACTORS





Retractors


Wundhaken
Separadores
Ecarteurs
Divaricatori



JOSEPH
BT-210/BT-214
160 mm




GILLIES
BT-15
180 mm



GILLIES
BT-16
180 mm



MINI-LANGENBECK
BT-186/BT-189
160 mm




Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori


Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori




Feine wundhaken flexibel-spitz
Delicate retractors flexible-sharp
Divaricatori delicati flessibili-acuti
Separadores finos flexibles-agudos
Écarteurs fins pointu flexible

Feine wundhaken flexibel-stumpf
Delicate retractors flexible-blunt
Divaricatori delicati flessibili-smussi
Separadores finos flexibles-romos
Écarteurs fins mousse flexible




Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori



DESMARRES
BT-70/BT-77
140 mm




BT-18

Spitz
Sharp
Acuto
Agudo
Pointu

BT-19

Stumpf
Blunt
Smusso
Romo
Mousse


SEÁN-MILLER
BT-18/BT-19
160 mm



BT-117
19 x 12 mm

BT-118
21 x 14 mm

KOCHER
BT-117/BT-118
205 mm




6,0 mm

10,0 mm



HURD
BT-1/BT-2
220 mm



12 mm

12 mm

Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori


Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori


Retractors
Divaricatori
Separadores
Ecarteurs




SENN GREEN
BT-48/BT-49
160 mm




RAGNELL
BT-27
140 mm



RAGNELL-DAVIS
BT-26
140 mm




KOCHER
BT-139
205 mm



KOCHER
BT-140
220 mm




LAHEY
BT-109
195 mm




CUSHING
BT-29
200 mm

Retractors


Wundhaken
Separadores
Ecarteurs
Divaricatori




BT-265



BT-264




BT-134
40 x 12 mm



BT-135
55 x 14 mm

DOCKHORN
BT-134/BT-135
210 mm




BT-264/BT-265
235 mm

Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori



Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori





Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori



OBWEGESER
BT-152/BT-157
220 mm





OBWEGESER
BT-92/BT-100
220 mm

Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori

Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori



Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori




Rostfrei - Stainless Steel
Acciaio - Acero - Acer


Biegsam - Malleable
Malleabili - Maleables
Souples

OM-85


Komplettset mit 3 Spatel
Set complete with 3 spatulas
Set completo con 3 spatole
Set completos con 3 espátulas
Set complete avec 3 spatules
OM-80 + OM-81 + OM-82




OM-80
200 mm



OM-81
200 mm




OM-82
240 mm



OM-87
200 mm



OM-88
200 mm



OM-100
180 mm

Schwarz Verchromt
Black anodized
Cromato nero
Cromado en negro
Chromé noir


Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori

Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori

Zungendrucker
Tongue Depressors
Abbassalingua
Bajalengua
Abaisse langue



Kinnhalter
Chin Holder
Divaricatore per mento
Separador para mentón
Écarteur pour menton





Fig. 1

Fig. 2


DO-08
140 mm

DO-09
150 mm

WIEDER
DO-08/DO-09



TOBOLD
OM-08
145 mm



HARTMANN
OM-06
150 mm




OBWEGESER
DO-40
155 mm

Retractors


Wundhaken
Separadores
Ecarteurs
Divaricatori



Ramushaken für Spina Nasalis
Ramus Retractor for nasal spine
Divaricatori per spina nasalis
Separadores para spina nasalis
Écarteurs pour spina nasalis



OBWEGESER
BT-67
220 mm




OBWEGESER
BT-68
220 mm

Mouth gags

Mundsperrer
Apribocca
Abrebocas
Écarteur buccal

Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori




ROSER KÖNIG
OM-60
185 mm




Paar Einsätzen aus Silikon für:
Pair of silicone inserts for:
Coppia di inserti in silicone per:
Pareja de inserto de silicone para:
Paire de inserts en silicon pour:

OM-40 - OM-41




OM-40R
30 mm

OM-41R
30 mm




HEISTER
OM-62
135 mm

Für Erwachsene - For adults
Per adulti - Para adultos
Pour adultes




MOLT
DOYEN-JANSEN
OM-41
140 mm

Für Erwachsene - For adults
Per adulti - Para adultos
Pour adultes



MOLT
DOYEN-COLLIN
OM-40
120 mm

Für kinder
For children
Per bambini
Para niños
Pour enfants



MOLT
OM-42
100 mm

Mouth gags

Mundsperrer
Apribocca
Abrebocas
Écarteur buccal




Abhalter für gesamten Ober-bzw. Unterkiefer
Retractor for complete upper resp. Lower jaws
Retratore per arcata superiore e inferiore intera
Separadores para quijadas enteras inferiores y superiores
Écarteurs pour machoires complètes inférieures et supérieures

OM-74 Fig. 0 100 x 95 mm


OM-75 Fig. 1 115 x 105 mm

OM-76 Fig. 2 115 x 110 mm


OM-77 Fig. 3 125 x 115 mm



ORINGER
OM-74/OM-77



MATY
OM-97
170 mm





OM-90
150 mm

Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori

Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori




Cheek Retractors

Wangenhalter


Divaricatori per guance

Separador para mejillas


Écarteurs des joues



STERNBERG
OM-50
160 mm




STERNBERG
OM-52
160 mm



COLUMBIA
OM-67
140 mm

Schwarz Verchromt
Black anodized
Cromato nero
Cromado en negro
Chromé noir


**NEW
BLACK**




BISHOP
OM-61
130 mm




MINNESOTA
OM-55
140 mm



CAWOOD-MINNESOTA
OM-66
160 mm



MINNESOTA
OM-55B
140 mm




CAWOOD-MINNESOTA
OM-66B
160 mm

Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori

Retractors


Wundhaken
Separadores
Ecarteurs
Divaricatori



Ramus Retraktoren Recht
Ramus Retractor Right
Retrattore Ramus destro
Separadores Ramus derecho
Écarteur Ramus droites


Ramus Retraktoren Link
Ramus Retractor Left
Retrattore Ramus sinistro
Separadores Ramus izquierdo
Écarteur Ramus gauches

8 mm



OBWEGESER
DO-41
165 mm

10 mm



OBWEGESER
DO-142
165 mm




MISCH
BT-158
210 mm



MISCH
BT-159
210 mm

Mouth gags

Mundsperrer
Apribocca
Abrebocas
Écarteur buccal




WHITEHEAD-JENNINGS

OM-136 90 mm

OM-137 110 mm

OM-138 130 mm

OM-139 150 mm



Komplett mit Rahmen, 3 Zungenspateln
und 2 beweglichen Seitenvälven
Complete with frame, 3 tongue depressors
and 2 adjustable lateral blades
Completo con telaio, 3 abbassalingua
e 2 valve laterali mobili
Completo con bastidor, 3 bajalenguas
y 2 valvas laterales ajustables
Complet avec cadre, 3 abaisse-langue
et 2 valves latérales réglables

DINGMANN


OM-110

Set completo

DINGMANN

Zungenspatel
Tongue depressor
Abbassalingua
Bajalengua
Abaisse langue


Fig. 1



OM-112

65 x 25 mm


Fig. 2



OM-113

70 x 29 mm

Fig. 3



OM-114


80 x 32 mm

Retractors


Wundhaken
Separadores
Ecarteurs
Divaricatori

Retractors


Wundhaken
Separadores
Ecarteurs
Divaricatori



GREENE
BT-131
225 mm



BT-132
200 mm



WASSMUND
BT-69
220 mm

Tongue Depressors - Tongue Holding Forceps

Zungenspatel - Zungenzangen


Abbassalingua - Pinze tiralingua

Bajalenguas - Pinzas tira-lengua

Abaisse langues - Pincers tire-langue




16 mm




OM-07
145 mm

19 mm




KROENIG OM-16
150 mm

18 mm




BUCHWALD OM-18
180 mm


23 mm



BRUENINGS OM-04
190 mm




FRENZEL OM-19
165 mm




OM-44


- 1 Paar Ersatzgummi
- 1 Pair of replacement rubber pads
- 1 Paio de gommini di ricambioi
- 1 Par de gomas de recambio
- 1 Paire de caoutchouc de rechange




YOUNG OM-45
160 mm



COLLIN OM-10/OM-30



OM-30
165 mm




OM-10
170 mm

Chain and Mouth Props

Kette und Mundstutzen
Catene e Apribocca
Cadena y Abrebocas
Chaîne et Ouvre Bouche

Retractors

Wundhaken
Separadores
Ecarteurs
Divaricatori



Dampfsterilisation bis 134°C
Autoclavable to 134°C
Autoclavabile a 134°C
Autoclave hasta 134°C
Autoclavable jusqu'à 134°C

Packung mit 2 Stück
Package with 2 units
Confezione da 2 pezzi
Caja con 2 unidades
Boîte avec 2 pièces

Schwarzer Gummi
Black Rubber
Gomma nera
Caucho negro
Caoutchouc noir



OM-70
SMALL
Klein, für Kinder
Small, for children
Piccolo, per bambini
Pequeño, para niños
Petit, pour enfants




OM-71
MEDIUM
Normal, für Erwachsene
Regular, for adults
Normale, per adulti
Regular, para adultos
Regulier, pour adultes



OM-72
LARGE
Groß, für Erwachsene
Large, for adults
Grande, per adulti
Grande, para adultos
Grand, pour adultes

IMPLANTOLOGY




Osteotomes

Osteotome
Osteotomi
Osteótomos
Ostéotomes


Implantology

Implantologie
Implantologie
Implantologia
Implantologiya



MINI-LAMBOTTE
FL-16/FL-58
125 mm

MINI-LAMBOTTE
FL-263/FL-269
170 mm



FREER
TG-10
160 mm

FREER
TG-11
160 mm

FREER
TG-12
160 mm

FREER
TG-20
160 mm


FREER
TG-21
160 mm

FREER
TG-22
160 mm


FREER
TG-23
160 mm

Osteotomes - Gouges

Osteotome - Meißel
Osteotomi - Sgorbie
Osteótomos - Gubias
Ostéotomes - Gouges




FL-115 6 mm
FL-114 8 mm
FL-116 10 mm




FL-132 6 mm
FL-135 8 mm
FL-136 10 mm

FL-114/FL-116
140 mm


FL-132/FL-136
140 mm




TG-15
170 mm




TG-30
155 mm




TG-31
155 mm



TG-32
155 mm



KRAMER-NEVINS
KN-55/KN-75
175 mm







DO-196
180 mm

Osteotomes - Gouges

Osteotome - Meißel
Osteotomi - Sgorbie
Osteótomos - Gubias
Ostéotomes - Gouges


Implantology

Implantologie
Implantologie
Implantologia
Implantologiya




Osteotomes - Chisels


Osteotome - Meissel
Osteotomi - Scalpelli
Osteótomos - Cinceles
Ostéotomes - Burins




MINI-LEXER
FL-59/FL-61
200 mm




MINI-LEXER
FL-230/FL-232
180 mm




TESSIER
FL-411 5 mm
FL-412 7 mm
160 mm




TESSIER
FL-413/FL-416
160 mm




EPKER
DO-180/DO-182
180 mm




EPKER
DO-183/DO-184
180 mm




EPKER
DO-185/DO-186
180 mm



COTTLE
DO-140
180 mm



COTTLE
DO-141
180 mm



COTTLE
DO-143
180 mm

Surgical mallets

Chirurgische Hämmer
Martelli chirurgici
Martillos quirúrgicos
Maillets chirurgicaux

Implantology

Implantologie
Implantologie
Implantologia
Implantologiya


Martelli chirurgici
Martillos quirúrgicos
Maillets chirurgicaux

Ø 20 mm - 140 g




FL-01
165 mm

Ø 25 mm - 300 g




**LUCAE
FL-10M**
190 mm

Ø 27 mm - 135 g



**HAJEK
FL-128**
200 mm


Ø 30 mm - 240 g



1 Seite flach - 1 Seite convex
1 Side flat - 1 Side convex
1 Lato piatto - 1 Lato convesso
1 Parte plana - 1 Parte convexe
1 Face plate - 1 Face convexe

**COTTLE
FL-104**
185 mm

Ø 20 mm - 160 g




FL-09R

Paar Ersatzbacken
Pair of replacement
Paioni ricambio
Par de recambio
Paire de rechange

**MEAD
FL-09**
190 mm

Mallet with interchangeable plastic disks
Hämmer mit auswechselbaren Kunststoffbacken
Maillet avec embouts de plastique interchangeables
Martillo con tapas plásticas reemplazables
Martello con punte di ricambi in teflon

Ø 25 mm - 200 g




FL-11R

Paar Ersatzbacken
Pair of replacement
Paioni ricambio
Par de recambio
Paire de rechange

**MEAD
FL-11**
200 mm

Mallet with interchangeable plastic disks
Hämmer mit auswechselbaren Kunststoffbacken
Maillet avec embouts de plastique interchangeables
Martillo con tapas plásticas reemplazables
Martello con punte di ricambi in teflon

Ø 26 mm - 320 g



FL-12R

Paar Ersatzbacken
Pair of replacement
Paioni ricambio
Par de recambio
Paire de rechange

**MEAD
FL-12**
190 mm


Raspatorien
Scollaperioste
Periostótomos
Rugines




Sinus Lift Retraktoren - Sinus Lift Retractors
Divaricatori per Sinus Lift - Separadores para Sinus Lift
Écarteurs pour Sinus Lift




SAUTTER
NK-249
145 mm




NK-105
125 mm




MODELL GIESSEN
NK-251
160 mm




ES-13
135 mm




FREER
180 mm



FREER
FR-16
185 mm



FREER-YASARGIL
FR-17
185 mm



OBWEGESER-FREER
718-215/718-216
220 mm

Trephines

Trepanbohrer

Trephine

Trépanos

Trépan

Implantology

Implantologie

Implantologie

Implantologia

Implantologiya

Combination Trocar-Trephine Instruments

Bestehend aus einem Trokar, einer gezahnten Trephine und einem Kolben zur Entnahme von Knochentransplantaten aus Tibia und Beckenkamm mittels eines geschlossenen Verfahrens

Consisting of trocar, a toothed trephine and a plunger for taking tibial and iliac bone grafts through a closed technique


Composto di un trequarti, una camicia dentata ed un pistone per il prelievo degli innesti ossei dalla tibia e dalla cresta iliaca tramite una tecnica chirurgica intracutanea

Compuesto por un trócar, un trépano dentado y un émbolo para obtener injertos óseos tibiales e iliacos mediante una técnica cerrada

Composé d'un trocart, d'une tréphine dentée et d'un piston pour le prélèvement de greffes osseuses de la tibia ou de la crête iliaque par une technique fermée

FK-100

Komplettset mit 4 Teilen
Set complete with 4 pieces
Set completo di 4 pezzi
Set completo con 4 piezas
Set complet avec 4 pièces



FK-105

Komplettset mit 4 Teilen
Set complete with 4 pieces
Set completo di 4 pezzi
Set completo con 4 piezas
Set complet avec 4 pièces

Implantology Instruments

Implantologie Instrumente

Strumenti per implantologia

Instrumentos para Implantología

Instruments pour Implantologie




MIR-30
130 mm



MIR-31
130 mm



MIR-32
130 mm




MIR-33
130 mm

BONE EXPANDERS

Distraktoren für Knochenleiste
Distracteurs de crête osseuse
Expansores de cresta ósea
Distrattori per cresta ossea

TUNNELING INSTRUMENTS


Tunnelierinstrumente | Instruments pour tunnellation | Instrumentos para tunelización | Strumenti per tunneling



UST-1/2ML



UST-1/2MS



UST-KPA



USS-YG7/8



UST-MP




UST-MPA

Osteotomes

Osteotome
Osteotomi
Osteótomos
Osteótomes

Implantology

Implantologie
Implantologie
Implantologia
Implantologiya



Osteotomes

Osteotome
Osteotomi
Osteótomos
Ostéotomes



# 1	Diam.
0	Ø 1.60
7.5	Ø 1.80
10	Ø 2.00
13.5	Ø 2.40
15.5	Ø 3.60

# 2	Diam.
0	Ø 1.90
7.5	Ø 2.20
10	Ø 2.35
13.5	Ø 3.30
15.5	Ø 3.40
23.5	Ø 4.50


# 3	Diam.
0	Ø 2.80
7.5	Ø 3.20
10.5	Ø 3.35
13.5	Ø 3.60
21	Ø 4.40

# 4	Diam.
0	Ø 3.10
7.5	Ø 3.70
10.5	Ø 3.90
15.5	Ø 4.30

# 5	Diam.
0	Ø 3.90
10	Ø 5.00
13	Ø 5.20
20.5	Ø 6.30

# 6	Diam.
0	Ø 5.50
10	Ø 6.20
13	Ø 6.30

1
Ø 1.6 mm




FL-120
155 mm

2
Ø 1.9 mm




FL-121
155 mm

3
Ø 2.8 mm



FL-122
155 mm

4
Ø 3.10 mm



FL-123
155 mm

5
Ø 3.90 mm



FL-124
155 mm

6
Ø 5.50 mm



FL-125
155 mm

Osteotomes

Osteotome
Osteotomi
Osteótomos
Ostéotomes


Implantology

Implantologie
Implantologie
Implantologia
Implantologiya

# 1	Diam.	# 2	Diam.	# 3	Diam.	# 4	Diam.	# 5	Diam.
0	Ø 1.50	0	Ø 2.10	0	Ø 2.60	0	Ø 3.00	0	Ø 3.50
5	Ø 1.90	5	Ø 2.60	5	Ø 3.00	5	Ø 3.50	5	Ø 3.90
8	Ø 2.20	8	Ø 2.80	8	Ø 3.30	8	Ø 3.80	8	Ø 4.20
10	Ø 2.40	10	Ø 3.00	10	Ø 3.50	10	Ø 4.00	10	Ø 4.40
12	Ø 2.60	12	Ø 3.20	12	Ø 3.70	12	Ø 4.20	12	Ø 4.60
14	Ø 2.80	14	Ø 3.40	14	Ø 3.90	14	Ø 4.40	14	Ø 4.80

Straight
Gerade
Droit
Recto
Diritto

Bone Pusher




Osteotomes

Osteotome
Osteotomi
Osteótomos
Ostéotomes




Meißeln Graduiert
Graduated Chisels
Scalpelli graduati
Cinceles graduados
Burins gradués

3.8 mm




FL-141
170 mm

3.8 mm




FL-142
170 mm

7.5 mm



FL-143
170 mm

7.5 mm




FL-144
170 mm

Osteotomes


Osteotome
Osteotomi
Osteótomos
Ostéotomes

Implantology


Implantologie
Implantologie
Implantologia
Implantologiya




FL-70/FL-73
165 mm



FL-74/FL-77
165 mm




FL-66/FL-69
165 mm



FL-65/FL-80
160 mm


Osteotomes

Osteotome
Osteotomi
Osteótomos
Ostéotomes



Dilatator
Dilator
Dilatatori
Dilatador
Dilataleur


A-Form



Spreader
Spreader
Expansori
Expansor
Expandeur

**D-Form
Convex Tip**

Ø 2 mm




FL-25
145 mm

Ø 3 mm



FL-26
145 mm

Ø 3.3 mm




FL-27
145 mm

Ø 4 mm



FL-28
145 mm

Ø 3 mm



FL-37
145 mm

Ø 4 mm



FL-38
145 mm

Ø 5 mm



FL-39
145 mm

Trephines


Trepanbohrer
Trephine
Trépanos
Trépans

Implantology

Implantologie
Implantologie
Implantologia
Implantologiya

For Tissue
Für Gewebe
Pour tissu
Para tejido
Per tessuto


Nicht Gezahnt
Non Serrated
No dentellato
No dentado
Non dentelé



Ø mm	GA-63	GA-64	GA-65	GA-66
Int.	2.5	4.0	4.5	5.5
Ext.	3.0	4.5	5.0	6.0


For Bone
Für Knochen
Pour os
Para hueso
Per osso

Gezahnt
Serrated
Dentellato
Dentado
Dentelé




Ø mm	GA-33	GA-34	GA-35	GA-36	GA-37
Int.	2.5	4.0	4.5	5.5	6.5
Ext.	3.0	4.5	5.0	6.0	7.0

Hautstanzen
Tissue Punches
Mucotomi Punches
Bisturi circular para tejido
Bistouri circulaire pour tissue




KEYES
FL-150/FL-153
100 mm



KEYES
FL-155/FL-159
90 mm

Implant Graduierung
Implant depth gauge with ball end
Sonda per controllare la profondità del sito
Sonda para controlar la profundidad



GA-60
100 mm


Implantology Instruments

Implantologie Instrumente

Strumenti per implantologia

Instrumentos para Implantología


Instruments pour Implantologie



**PELTON-CRANE
TS-40**
Connection Piece
int. Ø 7.5 mm
ext. Ø 10.0 mm
75 mm

inox Stainless Steel

Schalen
Bowl
Ciotola
Cuenco
Bacinella



- CS-01** Ø 40 mm, 15 cc - 0.5 oz
- CS-02** Ø 60 mm, 75 cc - 2.5 oz
- CS-03** Ø 80 mm, 225 cc - 7.5 oz
- CS-04** Ø 100 mm, 325 cc - 11 oz
- CS-05** Ø 120 mm, 500 cc - 17 oz
- CS-06** Ø 150 mm, 600 cc - 20 oz




Knochenschale
Bone Well
Bacinella per osso
Cubeta para hueso
Cuvette à os

Amalgambecher
Amalgam Well
Mortajo per amalgama
Recipiente para amalgama
Gobelet pour amalgame

SS-65
Ø 35x35 mm
depht 15 mm
5cc - 0.17oz

Absaugkanülen für Knochen mit Filter
Suction tube for bone with filter
Cannule d'aspiration avec trattenuta osso con filtro
Cánulas aspiradoras para hueso con filtro
Canule d'aspiration pour os avec filtre



Ersatzfilter
Filter of replacement
Filtro di ricambio
Filtro de recambio
Filtre de rechange

Ø 10 mm int.
Autoclavable
Stainless Steel



GF-134
20 mm

GF-133
180 mm

Membranplatzierungsinstrumente
Membrane Placement Instruments
Strumenti per posizionare membrana
Instrumentos para colocar membranas
Instruments pour application de membrane




MPI-1
185 mm




2/1



1/1



**LUNIATSCHEK
PF-6**
170 mm



1/1



1/1


Aufbring Instrument für Knochenmaterial
Applicator for bone material
Applicatore per materiale osseo
Cucharillas para material osoo
Applicateur pour matière osseuse



IM-15
180 mm



1/1



1/1

Bone Mills and Syringes

Knochenmühlen - und spritzen
Macina e siringhe per osso
Molinos y jeringas para hueso
Moulins et seringues à os

Implantology


Implantologie
Implantologie
Implantologia
Implantologiya




GA-18
Lower
Titanium




GA-17
Upper
Titanium




GA-19
Handle
with screw



GM-01
Bone Mill / Large
Dimensions: Ø 38 x 72 mm
Mill hole: Ø 0.71 mm
Volume: 2 cc.
Particle size: 710 µ (microns)




GM-02
Bone Mill / Small
Dimensions: Ø 30 x 44 mm
Mill hole: Ø 0.65 mm
Volume: 1 cc.
Particle size: 650 µ (microns)




GA-20
180 mm
Bone Crusher
Forceps

Ø int. 2.5 mm
Ø ext. 3.5 mm




SS-60M
160 mm

Ø int. 3.5 mm
Ø ext. 4.5 mm




SS-61M
160 mm

Ø int. 4.4 mm
Ø ext. 5.5 mm



SS-62M
160 mm

Ø int. 6 mm
Ø ext. 7 mm



SS-64M
160 mm


Bone Tampers - Bone Grafting Packers

Knochenstopfer

Compattatore per osso

Compactador para hueso

Compacteur à os




Bone-Grafter. Stainless Steel


Knochenschaber. Rostfrei Stahl
Grattini per osso. Inox
Raspador de hueso. Acero inoxidable
Racleur d'os. Acier inoxydable

Implantology

Implantologie
Implantologie
Implantologia
Implantologiya



Knochenschaber mit Reibeklinge
Bone-Grafter with blade
Grattino con lama
Raspador de hueso con hoja
Racleur d'os avec lame



Measuring instruments

Messinstrumente


Strumenti per misura

Instrumentos para medir

Instruments à mesurer




Meßbereich: 0 - 20 mm
Measuring range: 0 - 20 mm
Campo di misura: 0 - 20 mm
Alcance de medición: 0 - 20 mm
Etendue de mesure: 0 - 20 mm




OR-19
CASTROVIEJO
90 mm


Meßbereich: 0 - 80 mm
Measuring range: 0 - 80 mm
Campo di misura: 0 - 80 mm
Alcance de medición: 0 - 80 mm
Etendue de mesure: 0 - 80 mm




0 - 40 mm
Meßbereich / Measuring range / Campo di misura
Alcance de medición / Etendue de mesure



CASTROVIEJO
OR-11/OR-13
90 mm



CASTROVIEJO
OR-12/OR-14
170 mm



CASTROVIEJO
OR-18
170 mm



Osteometer

Osteometer
Osteometro
Osteometro
Osteometer

Implantology

Implantologie
Implantologie
Implantologia
Implantologiya

TRIDIMENSIONAL IMPLANTOLOGY CALIPER




DB-300


Schieblehre für Implantologie / Trei-Dimensional
Implantology Calliper / Tridimensional
Calibro per Implantologie / Tridimensionale
Calibre para Implantologia / Tridimensional
Calibre pour Implantologie / Tridimensionnel

BONE BLOCK CLAMP


Bone Holding Clamps / Knochenhalteklammern / Daviers à os
Pinzas de sujeción para hueso / Pinze per fissare l'osso



P-340
220 mm



P-341
200 mm



P-342
200 mm

Gewinkelt
Angled
Angolata
Curve
Courbé

Gewinkelt
Angled
Angolata
Curve
Courbé

Retractors - Bone Holding Forceps

Retraktoren - Knochenhaltezange


Divaricatori - Pinza per fissare l'osso

Separadores - Pinzas para sujeción para hueso


Écarteurs - Daviers à os




Tibia Wundspreizer
Tibial and Iliac Retractor
Divaricatore per tibia e cresta iliaca
Separadores tibiales e iliacos
Écarteur pour tibia e crête iliaque




**WEITLANER
BV-23**
165 mm



**RADOLF
FO-159**
135 mm




**WOLFF
LX-113**
180 mm



**DINGMANN
FO-130**

**DINGMANN
Mod. Fenestrated Jaw
FO-131**

**DINGMANN
FO-130/FO-131**
190 mm



**DINGMANN - Mod. Fenestrated Jaw
FO-131S**
190 mm

CURETTES & BONES





Periosteals (Raspatories)

Raspatorien

Scollaperiostei (Periostotomi)

Periostótomos (Legras)





Rugines périostiques



PERIOSTEAL ELEVATORS

PERIOSTEAL RASPATORIES

Scharf - Sharp - Taglienti - Cortante - Tranchant



OBWEGESER
DO-59
175 mm

OBWEGESER
DO-63
175 mm

OBWEGESER
DO-65
175 mm

OBWEGESER
DO-67
175 mm

WILLIGER
WG-10
165 mm

WILLIGER
WG-12
160 mm


WILLIGER
WG-14
140 mm

Periosteals (Raspatories)

Raspatorien
Scollaperiostei (Periostotomi)
Periostótomos (Legras)
Rugines périostiques

Periosteals / Curettes / Bone

Periost / Küretten / Knochen
Périosteales / Curettes / Os
Periostali / Curette / Osso
Nadkostnitsy / Kyuretki / Kost'



Raspatories - Elevators

Raspatorien - Elevatoren


Periostotomi - Elevatori

Legras - Elevadores

Rugines - Elevateurs




Rugines - Elevateurs




Halbscharf
Semi-Scharp
Semitagliente
Semicortante
Semitranchant

HENAHAN
BT-05/BT-07
175 mm

HENAHAN
BT-07M
175 mm



HENAHAN
BT-183/BT-185
200 mm



Raspatories - Elevators

Raspatorien - Elevatorien

Periostotomi - Elevatori

Legras - Elevadores

Rugines - Elevateurs


Periosteals / Curettes / Bone

Periost / Küretten / Knochen

Périosteales / Curettes / Os

Periostali / Curette / Osso

Nadkostnitsy / Kyuretki / Kost'




Raspatories - Elevators

Raspatorien - Elevatorien

Periostotomi - Elevatori

Legras - Elevadores

Rugines - Elevateurs



Raspatories - Elevators

Raspatorien - Elevatorien

Periostotomi - Elevatori

Legras - Elevadores

Rugines - Elevateurs


Periosteals / Curettes / Bone

Periost / Küretten / Knochen

Périosteales / Curettes / Os


Periostali / Curette / Osso

Nadkostnitsy / Kyuretki / Kost'



Legras - Elevadores
Rugines - Elevateurs

Instrumente für Sinuslift / Sinus Lift Instruments / Strumenti per rialzo del seno
Instrumentos para elevación de seno maxilar / Instruments d'élévation sinusienne




Raspatories - Elevators

Raspatorien - Elevatorien

Periostotomi - Elevatori


Legras - Elevadores

Rugines - Elevateurs



Instrumente für Sinuslift
Sinus Lift Instruments
Strumenti per rialzo del seno
Instrumentos para elevación de seno maxilar
Instruments d'élévation sinusienne

Fig. 24G



Raspatories - Elevators

Raspatorien - Elevatorien

Periostotomi - Elevatori

Legras - Elevadores

Rugines - Elevateurs


Periosteals / Curettes / Bone

Periost / Küretten / Knochen


Périosteales / Curettes / Os

Periostali / Curette / Osso


Nadkostnitsy / Kyuretki / Kost'




TASO-1



TASO-2




TASO-3




TASO-4




TASO-5



TASO-6



TASO-7



TASO-8


Raspatories

Raspatorien
Periostotomi
Legras
Rugines



Curettes
Küretten
Curettes
Curetas
Curettes

Instrumente für Sinuslift
Sinus Lift Instruments
Strumenti per rialzo del seno
Instrumentos para elevación de seno maxilar
Instruments d'élévation sinusienne




Curettes


Küretten
Curettes
Curetas
Curettes

Periosteals / Curettes / Bone

Periost / Küretten / Knochen
Périosteales / Curettes / Os
Periostali / Curette / Osso
Nadkostnitsy / Kyuretki / Kost'



Runder Hohlgriff 8 mm
Hollow Round Handle 8 mm
Manico tondo vuoto 8 mm
Mango hueco 8 mm
Manche rond vide 8 mm




720-101
180 mm

720-102
180 mm

720-103
180 mm

720-104
180 mm

HEMINGWAY
720-101/720-104
180 mm



720-85
85

720-86
86

720-87
87

720-88
88

720-85Z
85
Serrated

720-86Z
86
Serrated

LUCAS
720-85/720-88
180 mm




Curettes


Küretten
Curettes
Curetas
Curettes




HEMINGWAY
C-1
170 mm




HEMINGWAY
C-2
170 mm



HEMINGWAY
C-3
170 mm



HEMINGWAY
C-4
170 mm



SL-21
165 mm

Wundhachen Mandibular
Mandibular Nerve Displace
Per spostamento
nervo mandibolare
Para nero mandibular

1/1


1/1

Curettes

Küretten
Curettes
Curetas
Curettes


Periosteals / Curettes / Bone

Periost / Küretten / Knochen
Périosteales / Curettes / Os
Periostali / Curette / Osso
Nadkostnitsy / Kyuretki / Kost'




KERPEL
FK-136
FK-138
145 mm


KERPEL
FK-137
FK-139
145 mm




DL-50
155 mm




DL-51
150 mm





DL-52
150 mm




FK-83
160 mm





JANSEN
FK-84
160 mm



UFFENORDE
FK-188
170 mm



UFFENORDE
FK-189
170 mm



WILLIGER
FK-197/FK-198
170 mm

Fig. 00
3,4 mm

1/1

FK-197


Fig. 1
5,2 mm

1/1

FK-198

Curettes

Küretten
Curettes
Curetas
Curettes



Curettes

Küretten
Curettes
Curetas
Curettes

Periosteals / Curettes / Bone

Periost / Küretten / Knochen
Périosteales / Curettes / Os
Periostali / Curette / Osso
Nadkoznitsy / Kyuretki / Kost'

Curettes

VOLKMANN
FK-60/FK-67
170 mm




Fig.0000
2,8 mm
1/1
FK-60

Fig.000
3,6 mm
1/1
FK-61

Fig.00
4,4 mm
1/1
FK-62

,Fig.0
5,2 mm
1/1
FK-63


Fig.1
6,8 mm
1/1
FK-64

Fig.2
8,5 mm
1/1
FK-65


Fig.3
10 mm
1/1
FK-66

Fig.4
11,5 mm
1/1
FK-67

COLLEGE
CURETTE
MIR-10
160 mm



BRUNS
FK-110/FK-116
165 mm



3 mm
1/1
FK-110

4 mm
1/1
FK-111

5 mm
1/1
FK-112


6 mm
1/1
FK-113

8 mm
1/1
FK-114

10 mm
1/1
FK-115

12 mm
1/1
FK-116

SCHEDE
FK-120/FK-126
170 mm



2,5 mm
1/1
FK-120

3,4 mm
1/1
FK-121

4,2 mm
1/1
FK-122

5,2 mm
1/1
FK-123

6,3 mm
1/1
FK-124

7,4 mm
1/1
FK-125

8 mm
1/1
FK-126

Surgical Curettes


Chirurgische Küretten

Curette chirurgiche

Curetas quirúrgicas


Curettes chirurgicales

Fig. 4




**MOLT
MO-4L**
195 mm

Fig. 2/4



**MOLT
MO-10**
180 mm

**MOLT
MO-4RT**
145 mm



**MOLT
MO-4LT**
145 mm






Fig. 2



**MOLT
MO-2**
145 mm

Fig. 4



**MOLT
MO-4**
145 mm


Bone Files

Knochenfeilen
Lime per osso
Limas para hueso
Limes à os

Periosteals / Curettes / Bone


Periost / Küretten / Knochen
Périosteales / Curettes / Os
Periostali / Curette / Osso
Nadkostnitsy / Kyuretki / Kost'

Fig. 67




HOWARD
DO-1
180 mm

Fig. 45




MILLER
DO-2
185 mm

Fig. 52




MILLER
DO-3
180 mm

Fig. 12




HOWARD
DO-4
170 mm

Fig. 12A



MILLER
DO-4C
170 mm

Fig. 64



MILLER
DO-5
175 mm

Bone Files


Knochenfeilen
Lime per osso
Limas para hueso
Limes à os



Runder Hohlgriff 8 mm
Hollow Round Handle 8 mm
Manico tondo vuoto 8 mm
Mango hueco 8 mm
Manche rond vide 8 mm



Parodontal Feilen
Periodontal Files
Lime parodontali
Limas periodontales
Limes parodontales




Tissue Nippers


Gewebezängchen
Tronchesi per gengive e tessuti
Pinzas para cortar tejidos
Pincers coupante pour tissue

Periosteals / Curettes / Bone


Periost / Küretten / Knochen
Périosteales / Curettes / Os
Periostali / Curette / Osso
Nadkostnitsy / Kyuretki / Kost'



CHOEN
LX-60
105 mm



ILIC
LX-50
115 mm



LX-52
140 mm




Bone rongeurs

Hohlmeißelzangen

Pinze ossivore

Pinzas gubias

Pinces gouges



Bone rongeurs

Hohlmeißelzangen

Pinze ossivore

Pinzas gubias

Pinces gouges


Periosteals / Curettes / Bone

Periost / Küretten / Knochen


Périosteales / Curettes / Os

Periostali / Curette / Osso


Nadkostnitsy / Kyuretki / Kost'




BLUMENTHAL
FO-08/FO-19



MINI FRIEDMAN
FO-15
120 mm



FRIEDMANN
FO-17/FO-221
140 mm



MEAD
FO-21/FO-216
165 mm


Bone rongeurs

Hohlmeißelzangen


Pinze ossivore

Pinzas gubias


Pinces gouges



LEMPERT
FO-25/FO-43
190 mm



LEMPERT
FO-229/FO-230
155 mm



LUER
FO-18
155 mm

Bone rongeurs

Hohlmeißelzangen

Pinze ossivore

Pinzas gubias

Pinces gouges


Periosteals / Curettes / Bone

Periost / Küretten / Knochen


Périosteales / Curettes / Os

Periostali / Curette / Osso

Nadkostnitsy / Kyuretki / Kost'




RUSKIN
FO-138/FO-139
150 mm



FO-138


FO-139




KLEINERT-KUTZ
FO-156/FO-157
155 mm

FO-156


FO-157



FRIEDMANN
FO-16
160 mm



BEYER
FO-20
180 mm



ZAUFGAL-JANSEN
FO-34
180 mm

Bone rongeurs

Hohlmeißelzangen


Pinze ossivore

Pinzas gubias


Pinces gouges




Pinzas gubias
Pinces gouges




CLEVELAND
FO-231 140 mm
FO-232 170 mm



CLEVELAND
FO-233 150 mm
FO-234 170 mm



GARDNER
FO-26/FO-77
145 mm






BÖHLER
FO-60/FO-61
150 mm

Bone cutting forceps

Knochensplitterzangen
Pinze ossivore
Pinzas cortantes para huesos
Pinces coupantes à os

Periosteals / Curettes / Bone

Periost / Küretten / Knochen
Périosteales / Curettes / Os
Periostali / Curette / Osso
Nadkostnitsy / Kyuretki / Kost'



PERIODONTICS ENDODONTICS RESTORATION



Explorers - Sickle Probes

Zahnsonden

Esploratori - Sonde falce

Sondas exploradoras

Sondes exploratrices


Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration

Parodontie / Endodontie / Restauration

Parodontologia / Endodonzia / Restauro

Parodontologiya / Endodontiya / Restavratsiya



SO-1
Fig. 23

SO-2
Fig. 11

SO-3
Fig. 6

SO-4
Fig. 9

SO-5
Fig. 2

SO-6
Fig. 10


Explorers - Sickle Probes

Zahnsonden


Esploratori - Sonde falce

Sondas exploradoras

Sondes exploratrices




SO-7
Fig. 12



SO-8
Fig. 6A



SO-9
Fig. 8A



SO-10
Fig. 17A



SO-11
Fig. 3A



SO-12
Fig. 3

Periodontal Probes

Parodontometer
Sonde parodontali
Sondas parodontales
Sondes parodontales



Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya



Do not use Ultrasonic Cleaning
Verwenden Sie keine Ultraschall-Reinigung
Ne pas utiliser nettoyage à ultrasons
No usar limpieza ultrasónica
Non usare lavatrici ad ultrasuoni




Do not use metal brushes
Keine Drahtbürste verwenden
Ne pas utiliser de brosses métalliques
No utilizar cepillos metálicos
Non utilizzare spazzole metalliche





CP-11
SM-11
160 mm




CP-12
SM-12
160 mm





CP-12S
SM-13
160 mm



CP-15UNC
SM-15
160 mm



WHO
SM-21
160 mm



CP-10
SM-19
160 mm

Periodontal Probes

Parodontometer

Sonde parodontali

Sondas parodontales


Sondes parodontales




Double-ended Explorer
Doppelendiger Zahnsonden
Sonde exploratrice double
Sonda exploradora doble
Sonda falce doppia





11/12 ODU
SO-20
165 mm





CATTONI
UT3/4
SM-20
165 mm



NABERS
2QN
SM-22
175 mm



WILLIAMS
23 CP-11
EX-11
165 mm




WILLIAMS
23 CP-12
EX-12
165 mm

Periodontal Probes

Parodontometer
Sonde parodontali
Sondas parodontales
Sondes parodontales

Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya




Endodontic Instruments

Instrumente für Endodontie

Strumenti per Endodonzia

Instrumentos para Endodoncia

Instruments pour Endodontie




Endodontic Instruments

Instrumente für Endodontie
Strumenti per Endodonzia
Instrumentos para Endodoncia
Instruments pour Endodontie

Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya

Spreaders



SP-1
Fig.1S

SP-2
Fig.2S

SP-3
Fig.D11

SP-4
Fig.MA57

SP-5
Fig.0

SP-6
Fig.00

SP-7
Fig.3

Cotton Applicators

Watteträger


Porta cotone

Portaalgodón


Applicateur de coton




Tamponstopfer
Gauze Packers
Porta tamponi
Portatampones
Porte tampons



LUNIATSCHEK
PF-6
170 mm



UNIVERSAL
PF-40
140 mm



UNIVERSAL
PF-41
180 mm

Mirrors and handles


Mundspiegel und Griffe
Specchietti e manici
Espejos y mangos
Miroirs et manches

Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya

DOUBLE SIDED FRONT SURFACE RHODIUM MIRROR

Stainless steel - Rhodium coating



Reflecting surface is on **FRONT** of the glass.
Clear image. No distortion. No double/ghost image.

Rhodium coating for scratch resistance.


Enhanced indirect vision & illumination of dark areas.

Ergonomic & stream-lined design for less bulk.

+Visibility in any position: **2 opposite reflective faces**

(no need to turn the mirror to change from right to left).

Smoothly finished & retractor shaft design for patient comfort.



Simple
Stem

RH-14

4 - Ø 22mm

Cone
Socket


RH-15A

5 - Ø 24mm

Box of 6

FRONT SURFACE RHODIUM MIRRORS

Stainless steel - Rhodium coating




Reflecting surface is on **FRONT** of the glass.
Clear image. No distortion. No double/ghost image.

Rhodium coating for scratch resistance.

Enhanced indirect vision & illumination of dark areas.

Used for retraction of lips, cheeks and tongue.



Simple
Stem

RH-4

RH-5

Box of 12

4 - Ø 22mm

5 - Ø 24mm

Cone
Socket

RH-4A

RH-5A

PLANE MIRRORS

Stainless steel - Aluminium coating



Re ecting surface is on **REAR** of the glass.
Aluminium coating. Good image quality.
Flat mirror produce double/ghost image.

Low-cost.



Simple Stem

Box of 12

P-3

P-4

P-5

3 - Ø 18mm

4 - Ø 22mm

5 - Ø 24mm

Concave with enlarging effect



Simple Stem

Box of 12

G-3


G-4

G-5

3 - Ø 18mm

4 - Ø 22mm

5 - Ø 24mm



MA-1

Simple
Stem
120 mm
Solid

MA-2

Simple
Stem
125 mm
Solid

AL-MA5

Simple
Stem
120 mm
Hollow

MA-3

Simple
Stem
125 mm
Hollow

MA-7

Simple
Stem
135 mm
Hollow

MA-1A
Cone
Socket
120 mm
Solid

MA-2A
Cone
Socket
125 mm
Solid

MA-7A
Cone
Socket
120 mm
Hollow

inox StainlessSteel

inox StainlessSteel

Cr Chromium

Al Aluminium

inox StainlessSteel


inox StainlessSteel

Mirrors and handles

Mundspiegel und Griffe
Specchietti e manici
Espejos y mangos
Miroirs et manches



Spiegeldurchmesser - Diameter of mirrors
Diametro dello specchio - Diamètre des miroirs
Diámetro de los espejos




Inox



Kehlkopfspiegel mit
Sechskantgriff
Laryngeal mirrors
with hexagonal handle
Specchietto laringeo
con manico esagonale
Espejo laringeo
con mango hexagonal
Miroirs laryngiens
avec manche hexagonal

G-8 Fig. 00
G-10 Fig. 0
G-12 Fig. 1
G-14 Fig. 2
G-16 Fig. 3
G-18 Fig. 4
G-20 Fig. 5
G-22 Fig. 6
G-24 Fig. 7
G-26 Fig. 8
G-28 Fig. 9
G-30 Fig. 10


G-8/G-30
180 mm



Spiegelgriff
Mirrors handle
Manico per specchietti
Mango para espejos
Manche à miroir

Verchromt
Chromium
Cromato
Cromado
Chromé

MA-6
120 mm



Kehlkopfspiegel
Laryngeal mirrors
Specchietto laringeo
Espejos laringeos
Miroirs laryngiens

P-000 Fig. 000
P-00 Fig. 0
P-0 Fig. 0
P-01 Fig. 1
P-02 Fig. 2
P-03 Fig. 3
P-04 Fig. 4
P-05 Fig. 5
P-06 Fig. 6


P-000/P-06
120 mm

Dental Syringes

Dental-Spritzen
Siringhe dentali
Jeringas dentales
Seringues dentaires

Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya




Dosing Lever Syringe
for intraligamentary anesthesia.
Delivering precisely 0.2 ml of anesthetic.
No soft-tissue numbness.


SD-7
Imperial Tip
1.8 ml

SD-7M
Metric Thread
1.8 ml

SD-9NA
Imperial Tip
1.8 ml
(Sold only in North America)



SD-7R
Autoclavable plastic
protection sleeve
for SD-7




Pen-style Syringe
for intraligamentary
anesthesia

SD-20
Imperial Tip
1.8 ml

SD-20M
Metric Thread
1.8 ml


Aspirating Anesthetic Syringe



Aspirating
syringes feature a
blunt harpoon rod.


Lightweight and
Ergonomic.

High quality
stainless steel.



SD-8
Imperial Tip
1.8 ml
140mm


SD-8M
Metric Thread
1.8 ml
140mm



SD-8W
Imperial Tip
1.8 ml
140mm

SD-8WM
Metric Thread
1.8 ml
140mm


Aspirating and Self Aspirating Anesthetic Syringe



Self-aspirating
mechanism allows to
aspirate by means of
counter pressure.

Lightweight and
Ergonomic.


High quality
stainless steel.



SD-8NA
Imperial Tip
1.8 ml
140mm
Silicone covered
ring handle
(Sold only in
North America)

Silicone
covered
ring handle

Self Aspirating Anesthetic Syringe



SD-88
Imperial Tip
1.8 ml
140mm


SD-88M
Metric Thread
1.8 ml
140mm

Dental Syringes

Dental-Spritzen
Siringhe dentali
Jeringas dentales
Seringues dentaires



Zylinder-Ampullen Spritzen
Cartridge Syringes
Siringhe per tubofiale
Jeringas de ampollas
Seringues d'ampoules




Curettes

Küretten
Curette
Curetas
Curettes

Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya




Scalers and Curettes

Scaler und Küretten

Scaler e Curette

Raspadores y Curetas

Détartreurs et Curettes



USG-1/2ET
Gracey
Curette 1/2ET



USG-3/4ET
Gracey
Curette 3/4ET



USG-5/6ET
Gracey
Curette 5/6ET



USG-7/8ET
Gracey
Curette 7/8ET



USG-9/10ET
Gracey
Curette 9/10ET



USG-11/12ET
Gracey
Curette 11/12ET



USG-13/14ET
Gracey
Curette 13/14ET



USG-15/16ET
Gracey
Curette 15/16ET

Scalers and Curettes

Scaler und Küretten

Scaler e Curette

Raspadores y Curetas

Détartreurs et Curettes


Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration

Parodontie / Endodontie / Restauration

Parodontologia / Endodonzia / Restauro

Parodontologiya / Endodontiya / Restavratsiya



USG-1/2
Gracey
Curette 1/2



USG-3/4
Gracey
Curette 3/4



USG-5/6
Gracey
Curette 5/6



USG-7/8
Gracey
Curette 7/8




USG-9/10
Gracey
Curette 9/10



USG-11/12
Gracey
Curette 11/12



USG-13/14
Gracey
Curette 13/14



USG-1/2S
Gracey
Curette 1/2S



USG-3/4S
Gracey
Curette 3/4S



USG-5/6S
Gracey
Curette 5/6S



USG-7/8S
Gracey
Curette 7/8S



USG-9/10S
Gracey
Curette 9/10S



USG-11/12S
Gracey
Curette 11/12S



USG-13/14S
Gracey
Curette 13/14S




Scalers and Curettes

Scaler und Küretten

Scaler e Curette

Raspadores y Curetas

Détartreurs et Curettes



Scalers and Curettes

Scaler und Küretten
Scaler e Curette
Raspadores y Curetas
Détartreurs et Curettes

Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya



USS-NV0
Nevi
Scaler
Nv0



USS-NV11/12
Nevi
Scaler
NV11/12



USS-NV13/14
Nevi
Scaler
NV13/14



USS-T2/3
Taylor
Scaler 2/3



USS-204S
Scaler
204S



USS-H6/7
Scaler
H6/H7



USS-15/33
Towner-Jaquette
Scaler
U15/J33



USS-5/33
Jaquette
Scaler
H5/J33



USS-15/30
Towner-Jaquette
Scaler
U15/J30



USS-J14/15
Jaquette
Scaler
J14/15




Scalers and Curettes

Scaler und Küretten


Scaler e Curette

Raspadores y Curetas

Détarteurs et Curettes



T-USG-1/2
Gracey
Curette 1/2




T-USG-5/6
Gracey
Curette 5/6




T-USG-7/8
Gracey
Curette 7/8




T-USG-11/12
Gracey
Curette 11/12




T-USG-13/14
Gracey
Curette 13/14




T-USS-B5/6
Barnhart
Curette 5/6



T-USS-YG7/8
Younger-Good
Curette YG 7/8



T-USS-204S
Scaler 204S



T-USS-G6/7
Scaler G6/G7


⚠ Use a dedicated sharpening stone for Titanium Scalers & Curettes.

Curettes - Files

Curette - Lima
Curetasc - Lima
Curettes - Limes
Curettes - Limes

Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya





Scalers

Scaler
Ablatori di tartaro
Limpiadores
Détartreurs



Scalers




Amalgam gun for retrograde


Amalgampistole für retrograde
Siringhe portamalgame retrogradi
Portaamalgama para retrógradas
Porte-amalgame pour rétrogrades

Periodontics / Endodontics / Restoration


Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya




SS-51
130 mm




SS-55
130 mm




SS-54
130 mm




SS-56
130 mm



SS-58
130 mm



SS-50
130 mm



SS-52
130 mm

Amalgam Carriers

Amalgamträger
Porta amalgama
Porta-amalgamas
Porte-amalgames



inox Stainless Steel


Amalgambecher
Amalgam Well
Mortaio per amalgama
Vaso para amalgama
Gobelet pour amalgame




Knochenschale
Bone Well
Bacinella per osso
Cubeta para hueso
Cuvette à os

SS-65

Ø 35 x 35mm
depht 15 mm
5cc - 0.17oz




SS-15



SS-16


Aus Metall
Metal tip
In metallo
De metal
En métal



SS-41
170 mm

Autoclavable

Aus Kunststoff
Plastic tip
In plastica
De plàstico
En plastique



SS-42
160 mm

Composite Instruments

Kompositinstrumente
Strumenti per composito
Instrumentos para composites
Instruments pour composites

Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya



Fig. FE1



FELT
OT-1

Fig. FE2



FELT
OT-2

Fig. FE3



FELT
OT-3

Fig. FE4



FELT
OT-4

Fig. FE5



FELT
OT-5

Fig. FE6



FELT
OT-6

Fig. GO1



GOLDSTEIN
OT-7

Fig. GO2




GOLDSTEIN
OT-8

Fig. GO3



GOLDSTEIN
OT-9

Fig. GO4



GOLDSTEIN
OT-10

Composite Instruments

Kompositinstrumente

Strumenti per composito

Instrumentos para composites

Instruments pour composites



Titanium-coated
Titan-beschichtet
Recouvrement en Titane
Recubrimiento de Titanio
Ricoperto in Titanio

Fig. CP1
1,5 mm



OT-11

Fig. CP2
2,0 mm




OT-12

Fig. CP3
2,0 mm




OT-13

Fig. CP4
2,5 mm




OT-14

Fig. CP5



OT-15

Fig. CP6



OT-16

Composite Instruments

Kompositinstrumente
Strumenti per composito
Instrumentos para composites
Instruments pour composites


Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya



FELT

GOLDSTEIN



FE-1

FE-2

FE-3

FE-4

FE-5

FE-6

FE-4/5

GO-1

GO-2

GO-3

GO-4


Cement Spatulas

Zementspatel

Spatole per cemento

Espátulas para cemento


Spatules à ciment



SR
140 mm



S-00
180 mm




S-0
185 mm




S-1
180 mm



S-2
180 mm



S-3
180 mm



S-4
200 mm

Filling Instruments

Füllungsinstrumente

Strumenti per otturazione

Instrumentos para obturación

Instruments pour obturation


Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration

Parodontie / Endodontie / Restauration

Parodontologia / Endodonzia / Restauro

Parodontologiya / Endodontiya / Restavratsiya



Amalgam Carvers

Amalgam Instrumente

Strumenti per modellare l'amalgama


Instrumentos para modelar las obturaciones

Instruments à modeler l'amalgame



Instrument für Drahrückzug
Instrument for retractor wire
Strumento per filo retrattore
Instrumento para hilo separador
Instruments pour fil séparateur

Cleoids - Discoids



CV-11

USC-2
USC

F-2/3
FRAHM

CD-1
Fig. 1 /2

CD-2
Fig. 4/5

CD-3
Fig. W.A.C.D.


CD-4
Fig. 3/6

Filling Instruments

Füllungsinstrumente
Strumenti per otturazione
Instrumentos para obturación
Instruments pour obturation

Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya



Filling Instruments

Füllungsinstrumente

Strumenti per otturazione


Instrumentos para obturación

Instrumentes pour obturation



Gerieft
Serrated
Zigrinati
Estrados
Striée


Glatt
Smooth
Lisci
Lisos
Lissée




OZ-1
Fig.1

OZ-2
Fig.2

OZ-3
Fig.3




OL-1
Fig.1L




OL-2
Fig.2L



OL-3
Fig.3L



**WARD
OC-1**




**WARD
OC-2**

Filling Instruments

Füllungsinstrumente
Strumenti per otturazione
Instrumentos para obturación
Instruments pour obturation

Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya



OP-1
Fig.1

OP-2
Fig.2

OP-3
Fig.3

OP-4
Fig.4

OP-5
Fig.5

OA-1
Fig.1

OA-2
Fig.2

OA-3
Fig.3

Filling Instruments

Füllungsinstrumente

Strumenti per otturazione

Instrumentos para obturación

Instruments pour obturation




OX-1




OX-2




OX-3




WESTCOTT
OS-1



WESTCOTT
OS-2



OM-2



ON-1


Dental Restoration

Restaurative Zahnheilkunde
Odontoiatria conservativa
Odontología conservadora
Odontologie conservatrice

Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya

Füllungsinstrumente
Filling Instruments
Strumenti per otturazione
Instrumentos para obturación
Instruments pour obturation



Filling Instruments

Füllungsinstrumente

Strumenti per otturazione


Instrumentos para obturación

Instruments pour obturation




MARKLEY PLUGGERS


Amalgam Instrumente
Amalgam Instruments
Strumenti per amalgama
Instrumentos para amalgama
Instruments à amalgame




PM-1
Fig. 1M




PM-2
Fig. 2M




PM-3
Fig. 3M



PM-4
Fig. 4M



PM-5
Fig. 5M




PM-6
Fig. 6M

Enamel Instruments


Schmelzinstrumente
Strumenti da smalto
Instrumentos para esmalte
Instruments à émail

Periodontics / Endodontics / Restoration


Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya




BLACK
TS-1
Fig.26



BLACK
TS-2
Fig.27



BLACK
TS-3
Fig.28



BLACK
TS-4
Fig.29

Excavators

Excavatoren
Escavatori
Excavadores
Excavateurs



E-0
Fig. 0



E-1
Fig. 1



E-2
Fig. 2




E-3
Fig. 3



E-4
Fig. 4



EC-1
Fig. 1



EC-2
Fig. 2


Periodontal Chisels

Meißel für Parodontologie
Scalpelli parodontali
Cinceles periodontales
Burins parodontaux

Periodontics / Endodontics / Restoration


Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya

Fig. 1



FEDI
OF-1
175 mm

Fig. 2




FEDI
OF-2
175 mm

Fig. 3



FEDI
OF-3
170 mm

Fig. 4



OCHSENBEIN
OSB-4A
170 mm

OCHSENBEIN
OSB-4
170 mm

Periodontal Chisels

Meißel für Parodontologie
Scalpelli parodontali
Cinceles periodontales
Burins parodontaux



Bevel IN
Fig.1



OCHSENBEIN-1
OSB-1
165 mm

Bevel OUT
Fig.2



OCHSENBEIN-2
OSB-2
165 mm



TG CHISELS
718-113



OCHSENBEIN-3
718-103



OCHSENBEIN-4
718-104



RHODES BACK ACTION
718-107

Gingivectomy Knives

Parodontalmesser

Bisturi per gengivectomia

Bisturíes para la gengivectomía

Bistouris pour la gengivectomie


Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration


Parodontie / Endodontie / Restauration

Parodontologia / Endodonzia / Restauro


Parodontologiya / Endodontiya / Restavratsiya




KIRKLAND
KKD-1
Fig.15K




KIRKLAND
KKD-2
Fig.16K




KIRKLAND
KKD-4
Fig.15K - Fig.16K



ORBAN
OBN-1
Fig.1



ORBAN
OBN-2
Fig.2



ORBAN
OBN-4
Fig.1 - Fig.2


Double-ended Explorer

Doppelendiger Zahnsonden

Sonda falce doppia

Sonda exploradora doble

Sonde exploratrice double




Double-ended Periodontal Probes

Doppelendiger Parodontometer
Sondes parodontales doubles
Sonde parodontali doppie
Sondas parodontales doubles


Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya



Double-ended Periodontal Probes

Doppelendiger Parodontometer
 Sondes parodontales doubles
 Sonde parodontali doppie
 Sondas parodontales doubles



Spatulas - Filling Instruments

Spateln - Füllungsinstrumente

Spatole - Strumenti per otturazione

Espátulas - Instrumentos para obturación

Spatules - Instruments pour obturation

Periodontics / Endodontics / Restoration


Parodontologie / Endodontie / Restauration

Parodontie / Endodontie / Restauration


Parodontologia / Endodonzia / Restauro

Parodontologiya / Endodontiya / Restavratsiya


Heidemann Spateln
Heidemann Spatulas
Spatole Heidemann
Espátulas Heidemann
Spatules Heidemann



HEIDEMANN
711-301
Fig.1




HEIDEMANN
711-302
Fig.2




HEIDEMANN
711-303
Fig.3

Ø 8 mm
Hollow
Hohl
Creux
Hueco
Vuoto




IPC Fig. 1
715-311
Interproximal



IPC/OA Fig. 2
715-312
Interproximal Off-Angle

Filling Instruments

Füllungsinstrumente
Strumenti per otturazione
Instrumentos para obturación
Instruments pour obturation



Filling Instruments

Füllungsinstrumente

Strumenti per otturazione

Instrumentos para obturación

Instrumentes pour obturation


Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration

Parodontie / Endodontie / Restauration


Parodontologia / Endodonzia / Restauro

Parodontologiya / Endodontiya / Restavratsiya



Filling Instruments

Füllungsinstrumente
Strumenti per otturazione
Instrumentos para obturación
Instruments pour obturation



Ø 8 mm
Hollow
Hohl
Creux
Hueco
Vuoto

712-PFI
170 mm

HOLLENBACK



712-161H



712-162H



712-163H



712-164H

MARKLEY PLUGGERS



714-101



714-102

Filling Instruments

Füllungsinstrumente

Strumenti per otturazione

Instrumentos para obturación

Instrumentes pour obturation


Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration

Parodontie / Endodontie / Restauration


Parodontologia / Endodonzia / Restauro

Parodontologiya / Endodontiya / Restavratsiya



Filling Instruments

Füllungsinstrumente
Strumenti per otturazione
Instrumentos para obturación
Instruments pour obturation



Amalgam carvers

Amalgam Instrumente

Strumenti per modellare otturazioni

Instrumentos para modelar las obturaciones

Instruments à modeler les obturations


Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration

Parodontie / Endodontie / Restauration


Parodontologia / Endodonzia / Restauro

Parodontologiya / Endodontiya / Restavratsiya



Cavity preparation instruments

Schmelzinstrumente
Strumenti da smalto
Instrumentos para esmalte
Instruments à émail



Amalgam carvers

Amalgam Instrumente

Strumenti per modellare otturazioni

Instrumentos para modelar las obturaciones

Instruments à modeler les obturations

Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration

Parodontie / Endodontie / Restauration


Parodontologia / Endodonzia / Restauro

Parodontologiya / Endodontiya / Restavratsiya

Instrument für Drahrückzug
Instrument for Retractor wire
Strumento per filo retrattore
Instrumento para hilo retractor
Instruments pour fil rattracteur


Cleoids - Discoids

Dykal Applikator
Placement Instruments
Applicatori per sottofondo
Aplicadores de fondo de cavidades
Applicateurs de fond de cavité



Excavators

Excavatoren
Escavatori
Excavadores
Excavateurs



Gingivectomy knives

Parodontalmesser

Bisturi per gengivectomia

Bisturíes para la gingivectomía

Bistouris pour la gingivectomie


Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration

Parodontie / Endodontie / Restauration


Parodontologia / Endodonzia / Restauro

Parodontologiya / Endodontiya / Restavratsiya



Curettes

Küretten
Curette
Curetas
Curettes




Curettes

Küretten
Curette
Curetas
Curettes


Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya



Curettes

Küretten
Curette
Curetas
Curettes



GOLDMAN-FOX
718-271
Fig. 1



GOLDMAN-FOX
718-272
Fig. 2



GOLDMAN-FOX
718-273
Fig. 3



GOLDMAN-FOX
718-274
Fig. 4



GOLDMAN-FOX
718-275
Fig. 5



GOLDMAN-FOX
718-277
Fig. 21

Scalers and Curettes

Scaler und Küretten

Scaler e Curette

Raspadores y Curetas

Détartreurs et Curettes


Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration

Parodontie / Endodontie / Restauration

Parodontologia / Endodonzia / Restauro

Parodontologiya / Endodontiya / Restavratsiya




Scalers and Curettes

Scaler und Küretten

Scaler e Curette

Raspadores y Curetas

Détartreurs et Curettes



TAYLOR
718-TO3
Fig. 2/3

718-204
Fig. 204

718-204S
Fig. 204S

718-204SD
Fig. 204SD

718-CIO3
Fig. 2/3

CRANE-CAPLAN
718-CK6
Fig. 6

YOUNGER-GOOD
718-YGO8
Fig. 7/8

Scalers and Curettes

Scaler und Küretten

Scaler e Curette

Raspadores y Curetas

Détartreurs et Curettes

Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration

Parodontie / Endodontie / Restauration

Parodontologia / Endodonzia / Restauro

Parodontologiya / Endodontiya / Restavratsiya



RULE 23/24
718-RO4
Fig. 3 /4



I. UNIVERSITY
718-268
Fig. 17/18



TOWNER-JAQUETTE
719-111
Fig. U15/33

Ø 8 mm
Hollow
Hohl
Creux
Hueco
Vuoto



TOWNER-JAQUETTE
719-125
Fig. H5/33



TOWNER-JAQUETTE
719-126
Fig. H6/H7


Rubber-Dam Instruments

Kofferdam Instrumente

Strumenti per diga


Instrumentos para dique de goma

Instruments pour digue de caoutchouc



AINSWORTH
DF-02
160 mm

Lochzangen
Punch Forceps
Pinze foradiga
Perforadores
Pinces à perforer




DF-02M
Feder mit Schraube
Spring with screw
Molla con vite
Muelle con tornillo
Ressort avec vis



DF-02R
Lochscheibe
Punch Disk
Disco divisore
Disco agujereado
Disque perforé



DF-02P
Stift, mit Schraube
Pin, with screw
Perno con vite
Punzòn con perno
Broche



IVORY
DF-00
160 mm

Lochzangen
Punch Forceps
Pinze foradiga
Perforadores
Pinces à perforer




DF-00M
Spiralfeder
Spiral Spring
Molla a spirale
Muelle espiral
Ressort en spirale



DF-00L + DF-00V
Feder mit Schraube
Spring with screw
Molla con vite
Muelle con tornillo
Ressort avec vis



DF-00R
Lochscheibe
Punch Disk
Disco divisore
Disco agujereado
Disque perforé



DF-00P + DF-00N
Stift, mit Schraube
Pin, with screw
Perno con vite
Punzòn con perno
Broche avec vis

Rubber-Dam Clamp Forceps

Kofferdam-Klammerzangen

Pinze per uncini

Pinzas para grapas dique de goma

Pinces porte-crampons


Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration

Parodontie / Endodontie / Restauration

Parodontologia / Endodonzia / Restauro

Parodontologiya / Endodontiya / Restavratsiya



BREWER
DF-10
170 mm

PALMER
DF-15
170 mm

STOKE
DF-16
170 mm

DF-18
170 mm

Matrix Retainers and Matrix Bands

Matrizenspanner und Matrizenbänder

Portamatrici e matrici

Portamatrices y bandas matrices

Porte-matrices et bagues matrices



Molaren, 12 Stück
Molars, 12 pieces
Molari, 12 pezzi
Molares, 12 piezas
Molaires, 12 Pièces

Premolaren, 12 Stück
Premolars, 12 pieces
Premolari, 12 pezzi
Premolares, 12 piezas
Prémolaires, 12 Pièces



DF-40/1
Fig.1



DF-41/4
Fig.4



DF-40/2
Fig.2



DF-41/5
Fig.5




DF-40/3
Fig.3




DF-41/6
Fig.6


Kofferdam-Klammerzangen
Rubber-Dam Clamp Forceps
Pinze per uncini
Pinzas para grapas dique de goma
Pincers porte-crampons




IVORY
DF-12
170 mm



IVORY
DF-43
Fig. 14



IVORY
DF-40
Fig. 1



IVORY
DF-41
Fig. 12

Rubber-Dam Instruments

Kofferdam Instrumente
Strumenti per diga
Instrumentos para dique de goma
Instruments pour digue de caoutchouc

Periodontics / Endodontics / Restoration

Parodontologie / Endodontie / Restauration
Parodontie / Endodontie / Restauration
Parodontologia / Endodonzia / Restauro
Parodontologiya / Endodontiya / Restavratsiya

Kofferdam-Klammern
Rubber-Dam Clamps
Uncini per diga
Clamps para diques
Crampons pour digue



IVORY
DF-126
26



IVORY
DF-127
27



IVORY
DF-128
28



IVORY
DF-129
29



IVORY
DF-130
30



IVORY
DF-131
31




IVORY
DF-151
51




IVORY
DF-160
60

Kofferdam-Rahmen
Rubber-Dam Frames
Telaio per diga
Cuadro para diques
Cadre à digue



DF-30
80x75 mm
Kinder
Children
Bambini
Niños
Enfants



DF-32
95x115 mm
Erwachsene
Adults
Adulti
Adultos
Adultes

Matrix Retainers and Matrix Bands

Matrizenspanner und Matrizenbänder


Portamatrici e matrici

Portamatrices y bandas matrices

Porte-matrices et bagues matrices



Molaren, 12 Stück
Molars, 12 pieces
Molari, 12 pezzi
Molares, 12 piezas
Molaires, 12 Pièces



DF-03/7
#7



DF-03/8
#8



DF-03/9
#9



DF-03/10
#10

Premolaren, 12 Stück
Premolars, 12 pieces
Premolari, 12 pezzi
Premolares, 12 piezas
Prémolaires, 12 Pièces



DF-03/11
#11



DF-03/12
#12



DF-03/13
#13

TOFFLEMIRE
Universal



DF-03/1
#1



DF-03/2
#2



DF-03/3
#3



DF-01/5 5 mm x 1 mt

DF-01/6 6 mm x 1 mt


DF-01/7 7 mm x 1 mt



DF-05/5 5 mm x 5 mt

DF-05/6 6 mm x 5 mt

DF-05/7 7 mm x 5 mt



DF-06 5 mm

DF-08 6 mm

DF-09 7 mm

IVORY
DF-06/DF-09
#8

Paar - Pair
Coppia - Par
Paire



#1


Rechts
Right
Destra
Derecha
Droite



#2

Links
Left
Sinistra
Izquierda
Gauche

NYSTRÖM
DF-04
#1 + #2




TOFFLEMIRE
DF-07
Universal

CROWN BRIDGE




Crown Instruments

Kroneninstrumente
Pinze per corone
Alicates para coronas
Pincers pour couronnes




**KERLING
DP-38A**
145 mm

Kupferringabnehmzange
Copper Ring Removing Pliers
Pinza per togliere anelli di rame
Pinza para sacar anillos de cobre
Pince à enlever les bagues en cuivre




**FURRER
DP-38**
145 mm




DP-38R
Schraube
Replacement Screw
Vite di ricambio
Tornillos de recambio
Vis de rechange


Teleskopkronenzangen - Pliers for telescope crowns
Pinze per corone telescopiche - Alicates para coronas telescópicas
Pincers pour couronnes telescopes



**BÖHM
DP-28**
150 mm



**BÖHM
DP-29**
150 mm



DP-24
120 mm



DP-24C
Schlüssel
Key
Chiavetta
Llave
Clé



DP-24R
Diamantierte Einsätze, paar
Diamond inserts, pairs
Punte diamantate, coppia
Puntas diamantadas, par
Pointes diamantées, pair


Crown Instruments

Kronen Instrumente
Strumenti per corone
Instrumentos para coronas
Instruments pour couronnes

Crown and Bridge

Krone und Brücke
Couronne et pont
Corona e Ponte
Korona i most


MRT-C7
Complete Bridge Extractor
MRT-7 + C-71 + C-72 + C-73 + C-74




Removing pliers
for temporary
bridges and crowns




Tungsten-Carbide Inserts
Mit Hartmetalleinsätzen
Plaquettes en carbure de tungstène
Plaquetas de carburo de tungsteno
Inserti in carburo di tungsteno




MRT/C7



MRT-7



MRT-6
170 mm



PH-60
170 mm

Crown Removers

Kronenentferner

Leva corone

Removedor de coronas


Enlève-couronne



Enlève-couronne

MRT-C3

Complete with 3 Attachments
C-34 + C-35 + C-36




C-34

C-35

C-36

MRT-C4

Complete with 3 Attachments
C-41 + C-42 + C-43




C-41

C-42

C-43

MRT-C1

Complete with 3 Attachments
C-11 + C-12 + C-13



C-11


C-12

C-13



C-30

Usando la punta C-30
e invertendo la molla
si usa come spingicorone



MRT-C3



MRT-C4




MRT-C1

Crown Removers

Kronenentferner
Leva corone
Removedor de coronas
Enlève-couronne

Crown and Bridge

Krone und Brücke
Couronne et pont
Corona e Ponte
Korona i most



**UNIVERSAL
DL-60**
150 mm




DL-C
140 mm



DL-V
130 mm



DL-0
130 mm



DP-136
Crown Pliers
140mm

TEMPORARY CROWN PLACEMENT KIT | DP-136

1 Crown Pliers + 20 Replacement Tips + 10g Adhesive Powder

Crown Pliers: designed to remove upper and lower crowns and bridges. An adjustable stopper allows to adjust the grasping tips to the width of the crown and locking into place. It will prevent any unwanted pressures to the crown avoiding potential damages. The double spring back action along with the hook-lock mechanism allows for a smooth grasping motion.

Polymer Plastic Replacement Tips: provide a perfect and safely grip for temporary ceramic crowns without damaging metal, resin, or ceramic surfaces. A better grip of the crown can be achieved if the tips are coated with the supplied adhesive powder which is dampened by saliva.

Adhesive Powder: denture adhesive used as "glue" for recementing crowns or bridges that have come off.

DP-137 / Replacement Kit for DP-136: 20 Replacement Tips + 10g Adhesive Powder



ENDURING PERFORMANCE

Beautifully Crafted

THE PURSUIT OF EXCELLENCE

At MAST PAK Surgical Corporation, the mission is clear: we dedicate ourselves to outstanding craftsmanship through the pursuit of excellence and attention to fine detail. A century later, Hugo Friedman's unerring commitment to quality and innovation lives on at MAST PAK Surgical Corporation. As a result, MAST PAK Surgical Corporation instruments are indisputably the most skillfully created precision orthodontic instruments in the world.



Every instrument we create becomes, in essence, a critical extension of each practitioner's skill. To achieve this, we handcraft each instrument individually, assuring instruments of the highest quality, durability and performance available to the profession.

PRODUCTION OF SUPERIOR PRECISION

The process is demanding and complex — as many as 30 different operations may be needed to produce a single MAST PAK Surgical Corporation Orthodontic instrument.

IMMUNITY STEEL

MAST PAK Surgical Corporation uses only the highest-quality forged steel with an optimal blend of chromium and carbon. These elements, combined with exclusive MAST PAK Surgical Corporation processing, assure the creation of an outstanding instrument that will be resistant to corrosion given proper use and care.

ORBIT JOINT

Keeps instrument tips aligned, ensuring proper functionality and frustration-free use.

LASER ENGRAVING

Each instrument is laser engraved with the applicable wire sizes.

SATIN FINISH

All MAST PAK Surgical Corporation Orthodontic instruments are finished by hand to a satin luster; this non-glare surface keeps your instrument looking beautiful.

ERGONOMIC DESIGN

The size, shape, handle length and weight are balanced for optimum comfort and confidence.

ORTHO INSTRUMENTS

SLIM COLLECTION 1

CUTTERS 3

UTILITY PLIERS 5

WIRE FORMING PLIERS 8

HAND INSTRUMENTS 12

MATHIEUS, SCISSORS & HEMOSTATS 15

ORTHO ACCESSORIES 18

Slim Collection

Slim-Kollektion
Collection Mince
Collezione sottile
Tonkaya kolleksiya



SLIM COLLECTION

less means more
for your practice.

These eight instruments feature a 10% to 50% smaller working end than their original designs. Less instrumentation means more intraoral adaptation, maneuverability and patient comfort, especially for children or patients with smaller mouths.



Slim Collection

Slim-Kollektion
Collection Mince
Collezione sottile
Tonkaya kolleksiya



Orthodontics

Kieferorthopädie
Orthodontie
Ortodonzia
Ortodontiya



SLIM FLUSH CUT & HOLD DISTAL END CUTTER

678-113

50% reduced working end enhances visibility and improves efficiency by eliminating the need for a cotton roll or piece of gauze to catch the cut portion.

WIRE SIZE

$\leq .012" \times .021" \times .025"$ / .30mm - .53mm x .64mm

U.S. Pat. 7,182,595



SLIM BAND REMOVING PLIERS

678-503

With a 21% thinner and more tapered working end compared to our original 678-207 model, the Slim Band Removing Pliers allow for easy accessibility helping limit patient discomfort.



SLIM MICRO CUTTER

678-500

The slender tip allows for easier access when removing elastics. The face of the cutter has been "scooped out" and reduced by 13% as compared to our original 678-110 to enhance visibility.

WIRE SIZE

$\leq .008" \times .012"$



SLIM DIRECT BOND BRACKET TWEEZERS

678-504

The 50% thinner working end as compared to our 678-212 give the practitioner a confident and more visible grasp when bringing the bracket to the patient's mouth. While the thinner back end of the instrument allows for an easy manipulation of the bracket's position, by a simple flip of the instrument, it allows for one versus two instruments to be used in bracket placement and manipulation.



ULTRA SLIM WEINGART PLIERS

678-501

As an extension of your hands, the serrated tips allow for better gripping & accessibility into tight spaces that are typically not possible with hands. These pliers help guide the archwire through buccal tubes and into brackets with ease with a 20% reduced working end as compared to our 678-201 model.

WIRE SIZE

$\leq .030"$



SLIM BUCCAL TUBE BONDER TWEEZERS

678-505

The ultra slim working ends of these tweezers are ideal for proper and effective buccal tube bond placement by allowing the practitioner to hold the buccal tubes reversely when the buccal tubes are to be bonded.



SLIM HAMMERHEAD NiTi PLIERS

678-502

Reduce a multi-step process down to one step with these wire-forming pliers. Simply squeeze and you are done; no heat required. With a 10% thinner working end than our 678-327 model, these less intrusive pliers are designed for ultimate patient comfort. For NiTi wire only.

NiTi WIRE SIZE

$\leq .025"$



DISTAL CINCH BACK HAND INSTRUMENT

678-506

Designed to bend up to .021" x .025" stainless steel archwire, the slender working ends help improve access and enhance visibility in the posterior region. Use this instrument to cinch back posterior stainless steel wire distal to the buccal tube.

WIRE SIZE

$\leq .021" \times .025"$

CUTTERS

defined by precision.

With our diamond-honed cutting edges and orbit-formed joint, you will experience sharp and accurately-aligned cutting edges for every procedure, every time.

For each individual area of focus, we provide a specifically engineered cutter. To meet these stringent standards, all of our blades are hand-finished to maximize cutting performance for the appropriate wire size.



**UNIVERSAL CUT & HOLD
DISTAL END CUTTER**

678-101
When used intraorally it safely holds the cut portion to be discarded. Patient comfort is enhanced by the smooth, tapered tip.

WIRE SIZE
≤ .012" to .021" x .025" /
.30mm - .53mm x .64 mm



**UNIVERSAL CUT & HOLD DISTAL
END CUTTER WITH LONG HANDLE**

678-101L
Same as 678-101, with long handle.

WIRE SIZE
≤ .012" to .021" x .025" /
.30mm - .53mm x .64 mm



**FLUSH CUT DISTAL END
CUTTER, NO HOLD**

678-102
An exceptional universal cutter for all types of wires.

WIRE SIZE
≤ .021" x .025" / .53mm x .64mm



HARD WIRE CUTTER, 15°

678-103
15° offset tip enhances access in difficult areas. Designed to cut hard wires.

WIRE SIZE
≤ .022" x .028" / .56mm x .71mm



HARD WIRE CUTTER

678-104
Designed to cut hard wires.

WIRE SIZE
≤ .022" x .028" / .56mm x .71mm



MULTI-USE CUTTER

678-105
Designed to cut soft wire pins and ligature wires.

WIRE SIZE
≤ .015" / .38mm



**MINI PIN & LIGATURE
CUTTER**

678-106
Slim profile is a plus in difficult to access areas. Designed to cut soft wire pins and ligature wires.

WIRE SIZE
≤ .012" / .30mm



**MICRO-MINI PIN & LIGATURE
CUTTER**

678-107
Ultra slim profile is a plus in the most difficult to access areas. Designed to cut soft wire pins and ligature wires.

WIRE SIZE
≤ .012" / .30mm



**MICRO-MINI PIN & LIGATURE
CUTTER WITH LONG HANDLE**

678-107L
Same as 678-107, with long handle.

WIRE SIZE
≤ .012" / .30mm



**MICRO-MINI PIN & LIGATURE
CUTTER, 15°**

678-109
15° offset tip with ultra slim profile is a plus in the most difficult to access areas. Designed to cut soft wire pins and ligature wires.

WIRE SIZE
≤ .012" / .30mm



MICRO CUTTER

678-110
Ligature cutter with reduced tip dimension. Designed to cut soft wire pins and ligature wires.

WIRE SIZE
≤ .012" / .30mm



**FLUSH CUT & HOLD
DISTAL END CUTTER**

678-111
Cuts and holds wire distally flush to the buccal tube while reducing debondings. Designed for all types of wires.

WIRE SIZE
≤ .020" / .51mm
U.S. Pat. 7,182,595



**SLIM FLUSH CUT & HOLD
DISTAL END CUTTER**

678-113
Cuts and holds wire distally flush to the buccal tube while reducing debondings. Designed for all types of wires.

WIRE SIZE
≤ .012" to .021" x .025" /
.30mm - .53mm x .64 mm
U.S. Pat. 7,182,595

Utility Pliers

Allzweckzange
Pinces utilitaires
Pinze utili
Universal'nyye ploskogubtsy



UTILITY PLIERS

putting power and precision directly in your hands.

Strong, yet lightweight and perfectly balanced, the dimensions of these instruments are perfectly proportioned, they fit your grasp with unsurpassed ease. Designed to provide superior control with effortless versatility, MAST PAK Surgical Corporation's Utility Pliers assure comfort to both you and your patients.



Utility Pliers

Allzweckzange
Pinces utilitaires
Pinze utili
Universal'nyye ploskogubtsy



Orthodontics

Kieferorthopädie
Orthodontie
Ortodonzia
Ortodontiya



WEINGART PLIERS

678-201

Versatile utility pliers with serrations. Tapered 32° angular design is convenient for placement and removal of archwires.



SLIM WEINGART PLIERS

678-202

Slim, tapered 20° tips with serrations offer greater access and versatility.



WEINGART PLIERS WITH LONG HANDLE

678-201L

Same as 678-201, with long handle.



HOW PLIERS

678-203

How utility pliers have serrated tips. Useful for placement and removal of archwires, pins and other auxiliaries.



HOW PLIERS, OFFSET

678-204

Similar to the straight How pliers, except tips are angled at 40° for accessing lingual and other hard to reach areas.



POSTERIOR BAND REMOVING PLIERS, SHORT

678-208


The offset alignment of the tip and occlusal pad is well suited for band removal.



BAND SEATING PLIERS

678-205

Serrated tips grip brackets to assist in placing bands.



FORCE MODULE SEPARATING PLIERS

678-210

Used for placement of elastic separators.



ADHESIVE REMOVING PLIERS

678-206

Reversible, replaceable blade is used to remove the bulk of composite material after debonding attachments.



COON STYLE LIGATURE TYING PLIERS

678-211

Most popular pliers for tying stainless ligature wires.

WIRE SIZE
≤ .015" / .38mm

Utility Pliers

Allzweckzange
Pinces utilitaires
Pinze utili
Universal'nyye ploskogubtsy



POSTERIOR BAND REMOVING PLIERS, LONG 678-207

The alignment of the tip and occlusal pad is well suited for posterior band removal.



TIP BACK PLIERS 678-216

Designed to bend down the distal extension of an archwire protruding past the buccal tube. A unique contra-angle design allows for easy access and bending wire intraorally. Available in long handle for easy reach to the posterior buccal tube.

WIRE SIZE
<= .022" x .025" / .53mm x .64mm



BRACKET REMOVING PLIERS 678-219

Works well for steel, ceramic and plastic brackets of all types.



SLIM CROWN & BAND CONTOURING PLIERS 678-221MC

Same as 678-221, with slimmer beak.



ANGULATED BRACKET REMOVING PLIERS 678-220L

High angle design in long handle is effective for removal of all bracket styles from anterior and posterior teeth.



STEINER TYING PLIERS 678-222

Ideal for tying stainless steel ligature wire.

WIRE SIZE
<= .015" / .38mm



CROWN & BAND CONTOURING PLIERS 678-221

Designed for fitting of stainless steel crowns.



BAND CRIMPING PLIERS 678-225

Thin tips allow for easy access into molar and bicuspid bands.



SLIM CROWN & BAND CONTOURING PLIERS 678-221M

Same as 678-221, with slimmer beak.



BAND CRIMPING PLIERS WITH SPRING 678-225M

Thin tips allow for easy access into molar and bicuspid bands.

Wire Forming Pliers

Drahtformzange
Pinces à former les fils
Pinza per formare fili
Shchipty dlya formovki provoloki



Orthodontics

Kieferorthopädie
Orthodontie
Ortodonzia
Ortodontiya

WIRE FORMING PLIERS

wire forming is more than a skill, it is an art.

These instruments perform as natural extensions of your professional artistry. Our impeccably smooth finished tips allow you to bend the wire without risk of damage. Whether you are shaping stainless steel or nickel titanium, MAST PAK Surgical Corporation has a hand-finished, precision instrument to help bring your brilliant workmanship to life.



Wire Forming Pliers

Drahtformzange
Pinces à former les fils
Pinza per formare fili
Shchiptsy dlya formovki provoloki



LIGHT WIRE BIRD BEAK WITH GROOVE

678-305
A longer cone than a standard bird beak is useful for forming precise loops in light wire. The groove allows for locating and duplicating loops.

WIRE SIZE
<= .020" / .51mm



JARABAK PLIERS
678-301

Versatile light wire bending pliers. Grooves assist with precise bending and closing of loops. Serrated surface grips wires.

WIRE SIZE
<= .020" / .51mm



TWEED LOOP FORMING PLIERS

678-306
Loop forming pliers with .045", .060" and .075" diameters for creating bull and helical loops in wires.

WIRE SIZE
<= .022" x .025" / .53mm x .64mm



THREE JAW PLIERS
678-302

Useful for contouring and bending wire.

WIRE SIZE
<= .030" / .76mm



ARCH BENDING PLIERS
678-307

Designed for placing first, second and third order bends, .050" blade width.

WIRE SIZE
<= .022" x .025" / .53mm x .64mm



OMEGA LOOP FORMING PLIERS

678-303
Graduated cone is designed to form precise omega loops: .045", .060" and .075".

WIRE SIZE
<= .022" x .025" / .53mm x .64mm



TWEED PLIERS, SLIM
678-307S

Designed for placing first, second and third order bends, .025" blade width.

WIRE SIZE
<= .022" x .025" / .53mm x .64mm



BIRD BEAK PLIERS
678-304

Versatile loop forming pliers for round wires.

WIRE SIZE
<= .030" / .76mm



RECTANGULAR ARCH BENDING PLIERS

678-308
Designed for placing first, second and third order bends, .070" blade width.

WIRE SIZE
<= .022" x .025" / .53mm x .64mm

Wire Forming Pliers

Drahtformzange
Pinces à former les fils
Pinza per formare fili
Shchiptsy dlya formovki provoloki



Orthodontics

Kieferorthopädie
Orthodontie
Ortodonzia
Ortodontiya



GENERAL CONTOURING PLIERS

678-315
Tapered, serrated pliers for ease of contouring archwires.

WIRE SIZE
<= .021" x .025" / .53mm x .64mm



LINGUAL ARCH FORMING PLIERS

678-309
Designed to create double back and triple back bends in .030" and .036" lingual arch wires.

WIRE SIZE
<= .036" / .91mm



MARCOTTE LOOPING PLIERS

678-316
For precise forming of helical loops and T loops. No right angles. Prevents stress fracture of wires.

WIRE SIZE
<= .022" x .025" / .53mm x .64mm



HOLLOW CHOP CONTOURING PLIERS

678-310
Perfect for subtle contouring of arch wires.

WIRE SIZE
<= .030" / .76mm



HOLLOW CHOP PLIERS

678-317
Same as 678-310, with reduced radius.

WIRE SIZE
<= .030" / .76mm



O'BRIEN PLIERS
678-311

For precise forming of small helical loops and intraoral tip back bends.

WIRE SIZE
<= .020" / .51mm



LIGHT WIRE BIRD BEAK
678-318

A longer cone than a standard bird beak is useful for forming precise loops in light wire.

WIRE SIZE
<= .020" / .51mm



MINI THREE JAW PLIERS
678-312

Tips are rounded for ease of contouring light wire.

WIRE SIZE
<= .030" / .76mm



NANCE LOOP FORMING PLIERS

678-319
For precise forming of loops in four lengths: 3, 4, 5 and 6 mm. Groove at base aids in closing loops.

WIRE SIZE
<= .022" / .53mm

Wire Forming Pliers

Drahtformzange
Pinces à former les fils
Pinza per formare fili
Shchipty dlya formovki provoloki



ADAMS PLIERS

678-320

Excellent precise right angle bends and appliance application.

WIRE SIZE

<= .040" / 1.01mm



OLD STYLE LOOP PLIERS

678-324

Provides optimum accuracy and precision in forming and closing loops with one cylindrical three-step beak and one triangular beak.

WIRE SIZE

<= .022" x .025" / .53mm x .64mm



V-STOP PLIERS

678-321

Places V bend in archwires, including nickel titanium, to prevent archwire from traveling.

WIRE SIZE

<= .022" x .025" / .53mm x .64mm



BIRD BEAK WITH CUTTER

678-325

Smooth bending surfaces for a range of loop sizes while providing a sharp blade for cutting.

WIRE SIZE

<= .022" x .025" / .53mm x .64mm



NiTi THREE JAW PLIERS

678-322

Designed for forming and contouring of all archwire, especially for NiTi.

WIRE SIZE

<= .020" / .51mm



LONG TAPERED BIRD BEAK PLIERS

678-326

Bends a variety of smooth loops.

WIRE SIZE

<= .022" x .025" / .53mm x .64mm



OPTICAL PLIERS

678-323

Excellent for forming helical closing loops and arch adjustments.

WIRE SIZE

<= .022" x .025" / .53mm x .64mm



HAMMERHEAD NiTi TIE-BACK PLIERS

678-327

Multi-use NiTi bender, no heat.

WIRE SIZE

<= .025" / .63mm

Hand Instruments

Handinstrumente
Instruments à main
Strumenti manuali
Ruchnyye instrumenty



Orthodontics

Kieferorthopädie
Orthodontie
Ortodonzia
Ortodontiya

HAND INSTRUMENTS

icons for the industry.

All of our Orthodontic hand instruments are available in a variety of ergonomic handle sizes to accommodate your personal preference. Finely tuned to achieve accurate weight, balance and tip alignment, MAST PAK Surgical Corporation produces only the highest caliber tips and handles, in the manner of precision and excellence you demand.



Hand Instruments

Handinstrumente
Instruments à main
Strumenti manuali
Ruchnyye instrumenty



AOT BRACKET LIGATOR
678-903



SMALL PLUGGER/
PICKER
678-907



SCALER/BAND PUSHER
678-904



LIGATURE DIRECTOR/
BAND PUSHER
678-908



LIGATURE DIRECTOR/
SCALER
678-905



LARGE PLUGGER/
PICKER
678-909



DULL SCALER/BAND
PUSHER
678-906



DULL SCALER
678-910

Hand Instruments

Handinstrumente
Instruments à main
Strumenti manuali
Ruchnyye instrumenty



Orthodontics

Kieferorthopädie
Orthodontie
Ortodonzia
Ortodontiya



DULL SCALER
678-911



DIRECT BOND BRACKET
HOLDER
678-212



LIGATURE DIRECTOR/
BAND PUSHER
678-912



SLIM DIRECT BOND
BRACKET HOLDER
678-212M



DULL SCALER
678-913



BAND PUSHER
BP300
Aids in placement and
burnishing of band.



SCALER/PICKER
678-914



BAND PUSHER WITH
SHORT TIP
BP300S
Same as BP300 except the
working end is shorter.



PICKER/LIGATURE
DIRECTOR
678-916



BAND PUSHER/SCALER
BPS1
Versatile instrument for band
placement and cleanup
of band cement.

Hand Instruments

Handinstrumente
Instruments à main
Strumenti manuali
Ruchnyye instrumenty



BAND PUSHER/SCALER

BPS2

Same as BPS1 except the working end of the band pusher is shorter.



LIGATURE DIRECTOR/ TUCKER

LDT



FILLION EXPLORER/ LIGATURE DIRECTOR

FE/LD

Ideal for placing and removing elastic ligatures with lingual brackets.



PLUGGER/LIGATURE PICKER

PLPK



LIGATURE DIRECTOR

LD1



LINGUAL TORQUING KEY

TK018

Designed for proper wire placement and improved torque expression in an .018" slot lingual system.



LIGATURE DIRECTOR

LDG2



STRAIGHT .020"

APPLICATOR

TK020S



LIGATURE DIRECTOR/ PLUGGER

LDP



SINGLE END

LIGATURE DIRECTOR

WBI-M/R

Mathieus, Scissors & Hemostats

Mathieus, Scheren und Hämostatika
Mathieus, Ciseaux & Hémostats
Mathieus, Forbici ed emostatiche
Mat'yes, nozhnitsy i krovoostanavlivayushchiye zazhimy



Orthodontics

Kieferorthopädie
Orthodontie
Ortodonzia
Ortodontiya

MATHIEUS, SCISSORS & HEMOSTATS

forged from Immunity Steel.

Created with the same meticulous attention to detail as our pliers, MAST PAK Surgical Corporation hinged instruments are also forged from Immunity Steel, making them resistant to rust and corrosion.



Mathieus, Scissors & Hemostats

Mathieus, Scheren und Hämostatika

Mathieus, Ciseaux & Hémostats

Mathieus, Forbici ed emostatiche

Mat'yes, nozhnitsy i krovoostanavlivayushchiye zazhimi



Ortho Accessories

Ortho-Zubehör

Accessoires orthopédiques

Accessori orto

Orto akssessuary



Orthodontics

Kieferorthopädie

Orthodontie

Ortodonzia

Ortodontiya

ORTHO ACCESSORIES



HALSTED-MOSQUITO
HEMOSTAT
H4



CROWN & GOLD
SCISSORS CURVED
SCGC



WIRE-CUTTING SCISSORS
WCS
Serrated blade; notch for soft
wire-cutting.



CROWN & GOLD
SCISSORS STRAIGHT
SCGS



ORTHO MATHIEU,
NARROW TIP PLIERS
678-330

Slimmer tip makes this ideal
for gripping and placing elastic
auxillaries.



ORTHO MATHIEU,
WIDE TIP PLIERS
678-331

With wider, stronger tips.



ORTHO MATHIEU,
PERMA SHARP
678-332

Thin tips made with
Carbide inserts.



REPLACEMENT PARTS

Replacement blade for 678-206
Quantity: 1
678-001

Replacement blade for 678-303
Quantity: 1
678-002

Silicone Insert for 678-111
Quantity: 3
678-011

Silicone Insert for 678-113
Quantity: 3
678-013

1/4" Replacement pads for 678-207, 678-208
Quantity: 6
678-014

1/4" Replacement pad for 678-206
Quantity: 6
678-014A

1/8" Replacement set screw for 678-303
Quantity: 1
678-018

3/16" Replacement pad for 678-207, 678-208
Quantity: 6
678-036